

Pablo Estevez<sup>1</sup>  
 Sabrina Flor<sup>2</sup>  
 Oriana Boscolo<sup>1</sup>  
 Valeria Tripodi<sup>1,3</sup>  
 Silvia Lucangioli<sup>1,3</sup>

## Research Article

# Development and validation of a capillary electrophoresis method for determination of enantiomeric purity and related substances of esomeprazole in raw material and pellets

<sup>1</sup>Department of Pharmaceutical Technology

<sup>2</sup>Department of Analytical Chemistry and Physicochemistry, Faculty of Pharmacy and Biochemistry, University of Buenos Aires, Buenos Aires, Argentina

<sup>3</sup>National Council of Scientific and Technical Research (CONICET)

Received April 3, 2013

Revised November 8, 2013

Accepted November 11, 2013

A capillary electrophoresis method using CDs for quality control of esomeprazole (ESO) in terms of enantiomeric purity and related substances in raw material and pellets was developed. ESO is the *S*-enantiomer of omeprazole (OMZ). Several parameters were evaluated, including type and concentration of buffer and CD, concentration of additives and electrolyte pH. Resolution between the enantiomers of OMZ obtained for each parameter tested was calculated and the presence of the main related substance such as OMZ sulfone was carefully monitored. The optimized system consisted of 100 mM Tris-phosphate buffer pH 2.5 with 20 mM 2-hydroxypropyl- $\beta$ -CD, 1 mM sodium dithionite, temperature at 15°C, voltage at 28 kV, and UV detection at 301 nm. Once optimized, the electrophoretic system was validated according to ICH guidelines. The limits of detection and quantification for R-OMZ were 0.6  $\mu$ g/mL (0.06% w/w of ESO) and 2.0  $\mu$ g/mL (0.2% w/w of ESO), respectively. A mean concentration of R-OMZ <0.2% limit established by the United States Pharmacopeia (USP) was found in the raw material and six-pellet samples of ESO. No other impurities were found in the samples under these conditions. Therefore, the developed method was found to be appropriate not only for enantiomeric quality control of ESO but also for the analysis of ESO and the main related substance in raw material and pharmaceutical formulations as well as for stability indicating studies.

### Keywords:

chiral capillary electrophoresis / enantiomeric purity / esomeprazole / omeprazole / related substances  
 DOI 10.1002/elps.201300334


## 1 Introduction

Esomeprazole (5-methoxy-2-[(*R*)-[(4-methoxy-3,5-dimethylpyridin-2-yl)methanesulfonyl]-1*H*-1,3-benzodiazole] (ESO) (Fig. 1) is a single optical isomer proton-pump inhibitor used in the treatment of acid-related diseases. It is a potent inhibitor of gastric acid secretion and it does not undergo chiral inversion in vivo [1]. ESO has demonstrated a better pharmacological effect than the racemic product, omeprazole (OMZ), due to its higher systemic absorption [2], since (*R*)-OMZ is stereoselectively hydroxylated by cytochrome P450 CYP2C19 enzyme [3]. ESO is a hydrophobic compound with weak basic properties and it can be degraded unless it is protected against acidic conditions [4]. Therefore, ESO is administered in gastroresistant formulations such as pellets or capsules.

Few analytical methods have been developed for the enantiomeric resolution of OMZ, using HPLC and CE as well. The United States Pharmacopeia (USP) official method for magnesium ESO involves the use of a HPLC with a CHIRALPAK<sup>®</sup> AGP column with UV detection for the determination of *R*-enantiomer impurity in raw material [5]. Other HPLC methods have already been reported for the enantioselective analysis of OMZ in biological samples or in bulk drug. These methods are based on chiral stationary phases, using mainly polysaccharide or protein-based columns [6–9]. The main disadvantage of the HPLC methods is that they usually involve the use of expensive chiral columns to cover a reasonably wide application range, which have a relatively short lifetime. Moreover, chiral HPLC systems require large volumes of organic solvents. Nevertheless, HPLC methods using chiral stationary phases were the preferred technique for this purpose until the advent of capillary electrophoresis in the last decade [10, 11]. CE has been successfully employed in chiral analyses for numerous enantiomeric compounds [12–15]. The main advantages of CE for this kind of analysis are the extremely high efficiency, instrumentation simplicity, low sample and reagent consumption, and speed in method development and analysis, particularly of charged, polar and chiral compounds, using chiral selectors [16]. Moreover, chiral CE has proven to be more versatile and much

**Correspondence:** Professor Silvia E. Lucangioli Department of Pharmaceutical Technology, Faculty of Pharmacy and Biochemistry, University of Buenos Aires, Junin 956, 6to piso, Buenos Aires CP1113, Argentina  
**E-mail:** slucangi@ffyb.uba.ar

**Abbreviations:** ESO, esomeprazole; OMZ, omeprazole; R-OMZ, (*R*)-omeprazole; USP, United States Pharmacopeia



**Figure 1.** Chemical structures of ESO (A), R-OMZ (B), and OMZ sulfone (C). The sulfur atom is the stereogenic center.

less expensive than the chiral HPLC methods [17]. The most used chiral selectors are the CDs, a family of torus-shaped malto-oligosaccharides have found many applications in recent years because of their ability to form inclusion complexes with a large number of molecules. Only a few reports on CE chiral separation of OMZ are presented, in this sense, separation of enantiomers in capsules of racemic OMZ (enantiomer relation *R/S* 1:1) was achieved using a BGE consisting of 40 mM phosphate buffer adjusted to pH 2.2, 30 mM methyl- $\beta$ -CD, and 5 mM sodium disulfide to inhibit the oxidation of OMZ in acid medium. This method was mainly applied for the determination of OMZ racemate [18]. Another enantioselective analysis of racemic OMZ in tablets was carried out using 3% sulfated  $\beta$ -CD in 20 mM phosphate buffer pH 4.0, but without the use of any antioxidant. Also in this work, a comparison between chiral HPLC and chiral CE for the determination of OMZ enantiomers was made and the results presented show that both methods are suitable for the analysis of OMZ, but the HPLC method using a chiral column is more sensitive and achieves better resolution of OMZ enantiomers [17].

In the present work, we have developed a simple CE method for the separation and quantification of ESO, its *R*-impurity (R-OMZ), and its main related substances (OMZ sulfone, SFN) in esomeprazol raw material and pharmaceutical formulations such as gastroresistant pellets, using CDs as chiral selectors. The proposed method was optimized for different experimental parameters and validated according to ICH guidelines. To the best of our knowledge, there are no reports of any successful method described for the determination of enantiomeric purity and related substances of ESO in raw material and finished product.

## 2 Materials and methods

### 2.1 Chemicals and reagents

The racemic drug OMZ and ESO standards were supplied by INAME (National Institute of Drugs, Argentina). Enteric-coated ESO pellets were supplied by Laboratorios Richet Argentina. OMZ sulfone, a main related substance, was ob-

tained from Sigma (St. Louis, MO, USA). For the purpose of our study, several CDs were tested as chiral selectors for the CE system. Sulfated  $\beta$ -CD,  $\beta$ -CD, (2,3,6-tri-*O*-methyl)- $\beta$ -CD, and 2-hydroxypropyl- $\beta$ -CD (1.0 molar substitution, average molecular weight = 1540, H107) were supplied by Sigma. Sodium dithionite, used to prevent the degradation of OMZ into its main degradation products under acidic conditions, as well as sodium hydroxide, potassium hydroxide, sodium borate, (2-amino-2-hydroxymethyl) propan-1,3-diol (Tris), phosphoric acid, and methanol (HPLC grade) were purchased from Merck (Darmstadt, Germany). Water was purified in an EASY Pure<sup>TM</sup> RF equipment (Barnstead, Dubuque, IA, USA).

### 2.2 Standard and sample preparation

OMZ, R-OMZ, SFN, and ESO stock solutions were prepared in methanol to achieve a concentration of 1 mg/mL. Working standard solutions were prepared by diluting the corresponding stock solutions in 20 mM pH 9.2, sodium borate buffer according to the USP official method of OMZ [19, 20] to achieve a final concentration of 100  $\mu$ g/mL for ESO, 5  $\mu$ g/mL for R-OMZ and SFN, these solutions were used during the optimization of the system and sample analysis. OMZ working standard solution was prepared by diluting to a final concentration of 100  $\mu$ g/mL. ESO samples were prepared by dissolving 100 mg of the enteric coated pellets (theoretical concentration: 21.8% w/w of ESO) with 2 mL of 0.1N NaOH, aided by 5 min of sonication, and then completed to a final volume of 25 mL with 20 mM sodium borate buffer pH 9.2. The resultant solutions were centrifuged and analyzed immediately. The final pH sample was carefully selected to avoid OMZ degradation [4].

### 2.3 Stress conditions

Oxidation: 25 mL of a 3% v/v hydrogen peroxide solution was added to 25 mg of ESO accurately weighted. Acidic: 20 mg of ESO with 100 mL of 0.5 M hydrochloric acid solution were refluxed for 1 h. Alkaline: similar condition to the method

described under “acidic” using 0.5 M sodium hydroxide solution. Light: 1 mg of ESO solution was exposed to white light for 1 wk. All samples were diluted with 20 mM sodium borate buffer pH 9.2 to obtain a final concentration of 100 µg/mL of ESO.

## 2.4 Equipment and analytical conditions

The electrophoretic separation was performed on a P/ACE™ MDQ (Beckman Coulter) Karat V.8 software with diode array detection set at 301 nm. The separation was performed in an uncoated fused-silica capillary with 60 cm of length (effective length to the detector: 50 cm) and 75 µm of internal diameter. Standard and samples were introduced into the capillary by applying 0.5 psi of pressure for 5 s (sample volume injected approximately 19.6 nL calculated by Beckman Coulter CE Expert Lite software [https://www.beckmancoulter.com]). Capillary and sample temperature were set at 15°C, each run was performed at normal polarity (28 kV). New capillaries were conditioned by rinsing with 0.5 M potassium hydroxide for 10 min. Before each analysis, the capillary was washed with 0.1 M potassium hydroxide (2 min), water (2 min), and the running buffer (2 min). The electrophoretic separations were carried out in 100 mM Tris-phosphate buffer pH 2.5, containing 20 mM 2-HPβCD as chiral selector and 1 mM sodium dithionite to prevent degradation of the enantiomers of OMZ during the analysis.

## 3 Results and discussion

### 3.1 Method optimization

The BGE was optimized to obtain the highest resolution with no degradation products under the conditions established. Type and concentration of buffer and CD, concentration of antioxidant and electrolyte pH were tested. Resolution between the enantiomers of OMZ obtained for each system tested was calculated and the presence of related substances was carefully monitored.

#### 3.1.1 Choice of chiral selector and concentration

Based on literature, various CDs were tested (β-CD, sulfated β-CD, trimethyl-β-CD, and 2-HPβCD) in a concentration range of 10–40 mM. The selection of a suitable CD for chiral separation was made according to the shape and size of OMZ molecule and the respective CDs. From all CDs tested, only 2-HPβCD was able to separate OMZ enantiomers. The effect of different concentrations of 2-HPβCD (10–40 mM) on the resolution and migration times was investigated. Chiral separation was achieved for all cases, but 20 mM was the concentration chosen as a compromise between resolution and speed of analysis. From 10 to 30 mM, migration times obtained for OMZ were very similar and 20 mM being the

**Table 1.** Resolution values between ESO and its *R*-enantiomer obtained during method optimization by evaluating 2-HPβCD concentration (pH value fixed at 2.5) and evaluating the influence of the pH of the BGE (2-HPβCD concentration constant at 20 mM)


2-HPβCD concentration (mM)	R-OMZ migration time (min)	ESO migration time (min)	Resolution
10	14.4	14.6	1.72
20	16.4	16.8	2.32
30	18.5	18.9	2.44
40	24.3	25.0	6.45
BGE pH			
7	7.1	7.1	0
6	8.2	8.2	0
5	9.4	9.4	0
4	14.2	14.2	0
3	16.1	16.4	2.15
2.5	16.4	16.8	2.32
2	>60	>60	ND

Resolution was calculated according to USP [25]. Values “0” for resolution mean no separation at all. ND stands for not determined.

concentration that gave the adequate resolution in a shorter time. Concentrations higher than 30 mM were difficult to achieve because of dissolution problems of the CD in the BGE. Therefore, the concentration of 30 mM was not chosen although it presented a minor improvement in resolution. At 40 mM, resolution between OMZ enantiomers increased significantly along with migration times (Table 1).

#### 3.1.2 Effect of buffer pH

Knowing that OMZ suffers degradation at pH values <7.4 [4], different systems were initially tested in a pH range of 8–11. At alkaline pH values, the ESO, and R-OMZ are not charged ( $pK_{a1} = 4.0$   $pK_{a2} = 8.8$ , <http://www.pulsus.com/cddw2001/PI.htm>), so it was necessary to include in the electrolyte a charged pseudostationary phase such as SDS (plus a CD) to allow OMZ to run separately from the EOF. However, the high magnitude of EOF did not allow appreciable differentiation between the electrophoretic mobility of both enantiomers. The influence of the EOF in the resolution has been discussed a lot of times in previously published reports. It has been shown that for cations, the resolution decreases with increasing EOF. Despite the gain in efficiency, the lack of selectivity leads to a loss in resolution [21]. Analyzing the values of  $\mu_{ap}$  (apparent mobility) obtained at pH 2.5, the difference between ESO and R-OMZ is of 2.5%. This difference becomes zero when the EOF increases 1.7 times at pH 4.0. Therefore, it was decided to work at an acidic pH value where the EOF is almost zero and both enantiomers are charged and migrate according to their own electrophoretic mobility. The pH was evaluated in a range of 2–7, 2.5 being the pH at which the highest resolution



**Figure 2.** (A) Electropherogram 1 shows ESO subjected to acidic conditions. It can be seen SFN (d), ESO (c) and an unknown product (a). Electropherogram 2 shows the identification of the compounds, including R-OMZ (b) using the standards available under normal conditions. BGE conditions: 100 mM Tris-phosphate buffer pH 2.5 with 20 mM 2-HP $\beta$ CD and 1 mM sodium dithionite. (B) Influence of pH on enantiomeric separation during method optimization. Separation of OMZ enantiomers is achieved from pH 3 to lower values. BGE conditions: 100 mM Tris-phosphate buffer pH 7 to 2.5 with 20 mM 2-HP $\beta$ CD and 1 mM sodium dithionite.

was obtained with no significant increase in migration times (Table 1). In this range, pH values from 7 to 4 did not allow separation of the enantiomers, while pH values from 3 to 2.5 were sufficient to achieve separation (Fig. 2).

### 3.1.3 Choice and effect of buffer concentration

The effect of buffer type on the enantioseparation of OMZ was assessed by performing electrophoretic separation with the chosen chiral selector. Sodium phosphate buffer of pH 4.0 at a concentration of 20 mM was initially tested, but no separation was achieved under these conditions. Moreover, several products of degradation were observed due to the acidic conditions, for example the sulfone derivate. To prevent degradation, 20 mM sodium borate buffer pH 9.2 was tested. In this case, no separation was achieved, probably due

to the high EOF. Finally, 100 mM Tris/phosphate buffer pH 2.5 was tested and separation was achieved but with the presence of degradation products, so a suitable antioxidant was selected to prevent it, since it seemed more likely to achieve separation only at low pH values due to the protonation of the pyridinium group in OMZ ( $pK_a = 4.0$ )  $< pH 4$ . The dissociation equilibrium of the weakly acidic benzimidazole group ( $pK_a = 8.8$ ) is not affected in this pH range. Full chiral separation was achieved only at 100 mM of buffer concentration. Lower concentrations gave shorter migration times but with less resolution. The concentration established was chosen as a compromise between resolution and speed, being resolution the priority.

### 3.1.4 Effect of antioxidant addition


The optimized system consisted of 100 mM Tris-phosphate buffer pH 2.5 with 20 mM 2-HP $\beta$ CD. However, degradation products were observed under these conditions, so to prevent degradation of ESO at this pH value, different antioxidants were tested by adding them to the BGE, sodium dithionite being the most effective (without affecting baseline stability) at a minimum concentration of 1 mM.

### 3.1.5 Effect of instrumental parameters on the resolution

Temperature and voltage parameters were investigated to evaluate the impact of these instrumental variations on enantiomer resolution. Temperature was tested in a range from 15 to 30°C whereas voltage was tested from 15 to 30 kV. An increase in temperature affected the resolution negatively since it may cause a decrease in the association constant for the enantiomer-CD complex and therefore, a decrease in enantioselectivity. A desirable migration time decrease was achieved with the temperature increment, but it is not enough to justify the lower resolution obtained. Voltage increase provided good resolution of OMZ enantiomers and shorter migration times with good peak shape. When 30 kV was tested, the current increased to an undesirable value, so a slight reduction was made to avoid high current values. A temperature of 15°C and a voltage of 28 kV (current: 110  $\mu$ A) were selected for further analysis as a compromise between resolution and migration times.

### 3.2 Migration order

The migration order was established by analyzing a sample of OMZ standard (100  $\mu$ g/mL) spiked with ESO (100  $\mu$ g/mL). R-OMZ migrates first, which could be considered an advantage since peak tailing of ESO at sample concentration (1 mg/mL) might occur. In the analysis of enantiomeric purity of ESO samples, R-OMZ determination could be impaired if this enantiomer migrated in second place (Fig. 3).



**Figure 3.** (A) Electropherogram of ESO pellets sample. *R*-enantiomer standard (5  $\mu\text{g/mL}$ ) can be seen on the electropherogram at the upper right corner along with ESO standard (1000  $\mu\text{g/mL}$ ). The *R*-enantiomer concentration in the sample is lower than the 0.2% limit. (B) OMZ standard at a concentration of 100  $\mu\text{g/mL}$ . Migration order was established, *R*-OMZ migrating first, then ESO. BGE conditions (A) and (B) 100 mM Tris-phosphate buffer pH 2.5 with 20 mM 2-HP $\beta$ CD and 1 mM sodium dithionite.

### 3.3 Method validation for ESO

Once optimized, the electrophoretic optimal system consisting of 100 mM Tris-phosphate buffer pH 2.5 with 20 mM 2-HP $\beta$ CD and 1 mM sodium dithionite was validated evaluating the specificity, linearity, LOD, LOQ, precision, accuracy, and robustness of the method according to ICH guidelines [22]. All validation results are shown in Table 2.

The ability of the method to indicate chemical stability was studied by accelerated stress conditions (acid, alkaline, oxidation, and light) and specificity was also examined by comparing electropherograms of excipient blanks of the pharmaceutical formulation tested (pellets). Degradation was observed when ESO was exposed to acidic conditions. The presence of OMZ sulfone was detected in the electropherograms (Fig. 2). Oxidation conditions also lead to significant degradation. Exposure to light stress and alkaline conditions resulted in minimal degradation of ESO.

Also, to detect impurities on the electrophoretic peaks of the samples corresponding to the enantiomers, peak homogeneity was measured and evaluated, using different combined techniques available in a DAD system [23]. Absorbance at two different wavelengths was measured and spectra from several peak sections were compared. Both techniques demonstrated that the peak corresponding to ESO in the samples present a high level of purity, with no interferences of the excipients involved in the pellets samples. Pharmaceutical formulation excipients of pellets samples did not interfere with the assay.

Linearity for ESO was evaluated in a concentration range from 20 to 200  $\mu\text{g/mL}$  ( $n = 5$ ). In the case of *R*-OMZ, SFN were prepared in a concentration range of 2.0–6.0  $\mu\text{g/mL}$  ( $n = 5$ ) in the presence of 1 mg/mL of ESO. Each standard solution was analyzed in triplicate and plots of peak area versus concentrations were constructed for each compound.

**Table 2.** Parameters of validation of the method for *R*-OMZ, SFN, and ESO

Parameter	<i>R</i> -OMZ			SFN			ESO		
Linear range ( $\mu\text{g/mL}$ )	2.00–6.00			2.00–6.00			20.0–200.0		
$R^2$	0.9843			0.9955			0.9912		
LOD ( $\mu\text{g/mL}$ )	0.6			0.2			1.4		
LOD (% w/w of ESO)	0.06			0.02					
LOQ ( $\mu\text{g/mL}$ )	2.0			0.6			4.6		
LOQ (% w/w of ESO)	0.20			0.06					
Precision (%RSD)									
Intraday ( $n = 6$ )									
Migration time	0.3			0.7			0.5		
Peak area	1.9			1.5			0.9		
Interday ( $n = 18$ )									
Migration time	0.8			1.0			0.8		
Peak area	1.5			1.3			1.7		
Accuracy									
Spiked levels (%)	80	100	120	0.2	0.3	0.5	0.2	0.3	0.5
Recovery in pellets (%)	95.5	97.5	98.8	97.5	98.9	98.9	97.9	00.4	99.8

**Table 3.** Percentage of recovery of R-OMZ enantiomer for samples prepared at different ratios of enantiomer concentration

R/S enantiomer proportion	R-OMZ recovery (%)	RSD (%) <i>n</i> = 3
0.01/1	99.2	1.7
0.05/1	99.7	1.3
0.10/1	100.4	0.7
0.20/1	98.9	1.2
0.50/1	100.8	0.9
1.00/1	98.3	1.2

The LOD and LOQ were determined, based on S/N. A relation of 3:1 was used for estimating the LOD, whereas a 10:1 relation was used for the LOQ. Limit concentrations of each impurity were expressed in microgram per microliter as well as percentage weight in weight of ESO.

Precision was evaluated for intraday (*n* = 6) and interday (*n* = 18) assay of OMZ R-enantiomer and SFN (2 µg/mL of R-OMZ and SFN with 1 mg/mL of ESO) and ESO (50 µg/mL) after applying the sample preparation procedure, and it was expressed as RSD for retention times and areas.

Accuracy was calculated from recovery studies. Placebo samples prepared with all excipients contained in the pellets formulation were spiked with ESO at 80, 100, and 120% concentration levels. Preparations of each level were assayed by triplicate. Percentages of recovery values were obtained in the range of 97.9–100.4% (Table 2). For R-OMZ and SFN, placebo samples were spiked with a concentration equivalent to 0.2, 0.3, and 0.5% w/w respect to ESO (1 mg/mL). Preparations of each level were assayed by triplicate. Percentages of recovery values were obtained in the range of 95.5–98.8% and 97.5–98.9%, respectively.

To identify possible sources of error when the specified conditions are slightly changed, robustness studies were conducted by making small but deliberate changes in electrophoretic conditions. The method's capacity to remain unaffected was studied for ESO samples and standard solutions, and parameters selected for two level variations (upper and lower) were buffer pH (2.5 ± 0.2), 2-HPβCD concentration (20 ± 2 mM), temperature (15 ± 2°C), and voltage (28 ± 2 kV). The effects of variations in the electrophoretic parameters were evaluated statistically applying Student's *t*-test according to international guidelines [24, 25] using migration time, mean theoretical plates, resolution between enantiomers and ESO, R-OMZ, and SFN concentration in samples. The data obtained in these experiments did not show a significant difference of means (*p*-value > 0.05). Only buffer pH was statistically significant at the levels evaluated for migration time. pH values < 2.5 considerably increased the migration time, whereas higher pH values decreased it, but compromising the resolution. Therefore, buffer pH should be carefully controlled during analysis. All other parameters did not affect the separation and determination of enantiomers concentration significantly.

**Table 4.** Comparison between the methods already available

Method	Detection	System				Validation data			Resolution		Product
		Chiral selector	Capillary effective length	Buffer or mobile phase	Antioxidant	Linear range (µg/mL)	LOD (µg/mL)	LOQ (µg/mL)	Precision intraday (%RSD)	Migration Area time	
HPLC [17]	UV 302 nm	Chiralpack AD column		Hexane:ethanol (40:60 v/v)	No	25.0–150.0	0.0061	0.020	1.7	3.3	Tablets
HPLC [9]	UV 299 nm	Chiralpack IA column		MRBE-EA-EtOH-DEA (60:40:50.1 v/v/v/v)	0.1% DEA	0.5–25	0.17	0.51	0.27	Unknown	Tablets
CZE [18]	UV 301 nm	30 mM methyl-β-CD	21 cm	40 mM phosphate pH 2.2	5 mM sodium disulfide	2.0–40.0	0.31	1.0	<2	1.54	Capsules
CZE [17]	UV 202 nm	3% sulfated-β-CD	53 cm	20 mM Phosphate pH 4.0	No	25.0–150.0	1.3	4.3	1.8	1.5	Tablets
CZE <sup>a)</sup>	UV 301 nm	20 mM 2-HP-β-CD	50 cm	100 mM Tris-phosphate pH 2.5	1 mM sodium dithionite	2.00–6.00	0.6	2.0	2.02	2.32	Pellets

<sup>a)</sup> "ND" stands for not determined.

a) Method proposed in this research.

The evaluation of stability of both OMZ enantiomers standard solutions (ESO and R-OMZ standard) at different time intervals (up to 6 h) in the autosampler presented a RSD value of 1.9% for the content of R-OMZ and a RSD value of 3.7% for ESO. These results show that there is no significant variation in the content of each enantiomer standard from time of preparation to conclusion of the analyses, so these standard preparations were considered stable for 6 h in the autosampler.

### 3.4 Application: Analysis of ESO pharmaceutical formulations

The electrophoretic system performance on pharmaceutical samples was evaluated for the determination of R-OMZ impurity in ESO raw material as well as enteric-coated pellets (Fig. 3). Preparation of standards and samples is shown in the Section 2.2 described previously. Triplicate determinations of the content of the R-OMZ in the samples were performed. A mean concentration below the 0.2% limit established by the USP for R-OMZ was found in the raw material and six-pellet samples [5]. No other impurities were found in the samples under these conditions. Moreover, the capability of the method in the analysis of pharmaceutical samples containing different ratios of R-OMZ was evaluated. *R/S* enantiomer proportions were prepared at 0.01/1, 0.05/1, 0.1/1, 0.2/1, 0.5/1, and 1/1. These samples were prepared by triplicate and analyzed. Percentages of recovery for these samples are presented in Table 3.

A comparison between the methods already available in literature and the method we propose is shown in Table 4. The reported HPLC methods present better resolution, precision and lower LOD and LOQ than the CZE methods. However, they involve the use of expensive columns and higher solvent consumption, making these methods much less versatile than CZE methods. Moreover, determination of related substances is not described. USP official method involves the use of three different HPLC systems for assay, chromatographic purity, and enantiomeric purity, respectively, while our CZE system allows performing all three determinations at the same time, which is a great advantage [5]. The other CZE methods reported [17, 18] present lower resolution values for enantiomer separation and do not show separation of related substances (sulfone and other degradation products).

## 4 Concluding remarks

The method developed is simple and robust for quality control of enantiomeric purity of ESO raw material and finished product allowing the resolution of the enantiomers and the main related substances such as the OMZ sulfone encoded in USP. This method has been validated for the analysis of raw material and pellet samples without any excipient or degradation products interference. To the best of our knowledge, this is the first method using CE applied to the determination of

R-OMZ impurity and its main degradation products in ESO pellets samples. This method is applicable for quality control of ESO products in pharmaceutical industry and a stability indicating method, being a valuable alternative to chiral analyses that employ HPLC, which are much more expensive and less versatile compared to chiral capillary electrophoresis.

*The authors acknowledge the financial support for this work provided by University of Buenos Aires.*

*The authors have declared no conflict of interest.*

## 5 References

- [1] Mckeage, K., Blick, S. K. A., Croxtall, J. D., Lyseng-Williamson, K. A., Keating, G. M., *Drugs* 2008, **68**, 1571–1607.
- [2] Saccar C. L., *Expert Opin. Drug Metab. Toxicol.* 2009, **9**, 1113–1124.
- [3] Tybring, G., Bottiger, Y., Widen, J., Bertilsson, L., *Clin. Pharmacol. Ther.* 1997, **62**, 129–137.
- [4] Mathew, M., Das Gupta, V., Bailey, R., *Drug Develop. Ind. Pharm.* 1995, **21**, 965–971.
- [5] United States Pharmacopoeia 32 Revision. *United States Pharmacopoeial Convention*, Rockville, MD, USA 2009, p. 2298.
- [6] Balmer, K., Persson, B. A., Lagerstrom, P. O., *J. Chromatogr. A* 1994, **660**, 269–273.
- [7] Tanaka, M., Yamazaki, H., Hakusui, H., *Chirality* 1995, **7**, 612–615.
- [8] Cass, O.B., Degani, A. L. G., Cassiano, N., *J. Liq. Chromatogr. Rel. Technol.* 2000, **23**, 1029–1038.
- [9] Zanitti, L., Ferreti, R., Gallinella, B., La Torre, F., Sanna, M. L., Mosca, A., Cirilli, R., *J. Pharm. Biomed. Anal.* 2010, **52**, 665–671.
- [10] Nishi, H., *Electrophoresis* 1999, **20**, 3237–3258.
- [11] Maier, N. M., Franco, P., Lindner, W., *J. Chromatogr. A* 2001, **906**, 3–33.
- [12] Tonon, M. A., Bonato P. S., *Electrophoresis* 2012, **33**, 1606–1612
- [13] Aguiar, F. A., de Gaitani C. M., Borges K. B., *Electrophoresis* 2011, **32**, 2673–2682.
- [14] Fakhari, A. R., Nojavan, S., Haghgoo, S., Mohammadi, A., *Electrophoresis* 2008, **29**, 4583–4592.
- [15] Guan, J., Yan, F., Shi, S., Wang, S., *Electrophoresis* 2012, **33**, 1631–1636.
- [16] Blaschke, G., Chankvetadze, B., *J. Chromatogr. A* 2000, **875**, 3–25.
- [17] Bonato, P. S., Paias, F. O., *J. Braz. Chem. Soc.*, 2004, **15**, 318–323.
- [18] Berzas Nevado, J. J., Castañeda Peñalvo, G., Rodriguez Dorado, R. M., *Anal. Chim. Acta.* 2005, **533**, 127–133.

- [19] United States Pharmacopoeia 31 Revision, *United States Pharmacopoeial Convention*, Rockville, MD, USA 2009, p. 2851.
- [20] Flor, S., Tripodi, V., Scioscia, S., Revello, L., Lucangiolli, S., *J. Liq. Chromatogr. Relat. Technol.* 2010, *33*, 1666–1678.
- [21] Kenndler, E., *J. Microcolumn Sep.* 1998, *10*, 273–279.
- [22] Guidance Q2 (R1), *Validation of Analytical Procedures: Text and Methodology. International Conference on Harmonization. Final version 2005*.
- [23] Farbe, H., Lebris, A., Blanchin, M., *J. Chromatogr. A* 1995, *697*, 81–88.
- [24] Dejaegher B., Vander Heyden Y., *J. Chromatogr. A* 2007, *1158*, 138–157.
- [25] United States Pharmacopoeia 32 Revision, *United States Pharmacopoeial Convention*, Rockville, MD, USA 2009, p. 227.



## Author Query Form

---

**Journal**      **ELPS**

**Article**      **elps201300334**

---

Dear Author


*During the copy-editing of your paper, the following queries arose. Please respond to these by marking up your proofs with the necessary changes/additions. Please write your answers clearly on the query sheet if there is insufficient space on the page proofs. If returning the proof by fax do not write too close to the paper's edge. Please remember that illegible mark-ups may delay publication.*

Query No.	Description	Remarks
Q1	Author: Please provide institute and town/city name in affiliations 1 and 3.	
Q2	Author: Both the terms R-OMZ and (R)-OMZ have been provided in the text. Please check for correctness.	
Q3	Author: The term "United States Pharmacopeia" has been abbreviated as USP in text. Please check for correctness.	
Q4	Author: Please check all heading levels as typeset for correctness.	
Q3	Author: Please check abbreviation list as typeset for correctness.	
Q5	Author: The term SFN has not been expanded in the text. Please check for correctness.	
Q6	Author: Few edits have been made to enhance the readability of the text. Please check for correctness.	
Q7	Author: "pka" has been changed to "pKa." Please check for correctness.	
Q8	Author: Please define CZE.	
Q9	Author: Please provide conference location in reference 22.	


USING e-ANNOTATION TOOLS FOR ELECTRONIC PROOF CORRECTION

Required software to e-annotate PDFs: **Adobe Acrobat Professional or Adobe Reader (version 8.0 or above).** (Note that this document uses screenshots from **Adobe Reader X**)  
 The latest version of Acrobat Reader can be downloaded for free at: <http://get.adobe.com/reader/>


Once you have Acrobat Reader open on your computer, click on the **Comment** tab at the right of the toolbar:



This will open up a panel down the right side of the document. The majority of tools you will use for annotating your proof will be in the **Annotations** section, pictured opposite. We've picked out some of these tools below:



**1. Replace (Ins) Tool – for replacing text.**




Strikes a line through text and opens up a text box where replacement text can be entered.


**How to use it**

- Highlight a word or sentence.
- Click on the **Replace (Ins)** icon in the Annotations section.
- Type the replacement text into the blue box that appears.

Standard framework for the analysis of microeconomics. Nevertheless, it also led to the emergence of strategic behavior. The number of competitors is that the structure of the industry, which led to the emergence of strategic behavior. The main components of the industry level, are exogenous and endogenous. The important works on entry by Cournot (1838) and henceforth, we open the 'black b



**2. Strikethrough (Del) Tool – for deleting text.**




Strikes a red line through text that is to be deleted.

**How to use it**

- Highlight a word or sentence.
- Click on the **Strikethrough (Del)** icon in the Annotations section.

there is no room for extra profits as the number of firms is zero and the number of firms (n) values are not determined by Blanchard and ~~Kiyotaki~~ (1987), perfect competition in general equilibrium. The effects of aggregate demand and supply in the classical framework assuming monopoly are an exogenous number of firms

**3. Add note to text Tool – for highlighting a section to be changed to bold or italic.**




Highlights text in yellow and opens up a text box where comments can be entered.


**How to use it**

- Highlight the relevant section of text.
- Click on the **Add note to text** icon in the Annotations section.
- Type instruction on what should be changed regarding the text into the yellow box that appears.

dynamic responses of mark-ups are consistent with the VAR evidence



**4. Add sticky note Tool – for making notes at specific points in the text.**




Marks a point in the proof where a comment needs to be highlighted.

**How to use it**

- Click on the **Add sticky note** icon in the Annotations section.
- Click at the point in the proof where the comment should be inserted.
- Type the comment into the yellow box that appears.

and supply shocks. Most of the time, the number of firms is zero and the number of firms is not determined by the standard framework. The effects of aggregate demand and supply in the classical framework assuming monopoly are an exogenous number of firms



USING e-ANNOTATION TOOLS FOR ELECTRONIC PROOF CORRECTION


**5. Attach File Tool – for inserting large amounts of text or replacement figures.**



Inserts an icon linking to the attached file in the appropriate place in the text.

**How to use it**

- Click on the **Attach File** icon in the Annotations section.
- Click on the proof to where you'd like the attached file to be linked.
- Select the file to be attached from your computer or network.
- Select the colour and type of icon that will appear in the proof. Click OK.



**6. Add stamp Tool – for approving a proof if no corrections are required.**




Inserts a selected stamp onto an appropriate place in the proof.

**How to use it**

- Click on the **Add stamp** icon in the Annotations section.
- Select the stamp you want to use. (The **Approved** stamp is usually available directly in the menu that appears).
- Click on the proof where you'd like the stamp to appear. (Where a proof is to be approved as it is, this would normally be on the first page).

of the business cycle, starting with the  
 on perfect competition, constant ret  
 production. In this environment goods  
 extra...  
 he...  
 determined by the model. The New-Ky  
 otaki (1987), has introduced produc  
 general equilibrium models with nomin  
 and...




**7. Drawing Markups Tools – for drawing shapes, lines and freeform annotations on proofs and commenting on these marks.**


Allows shapes, lines and freeform annotations to be drawn on proofs and for comment to be made on these marks..

**How to use it**

- Click on one of the shapes in the **Drawing Markups** section.
- Click on the proof at the relevant point and draw the selected shape with the cursor.
- To add a comment to the drawn shape, move the cursor over the shape until an arrowhead appears.
- Double click on the shape and type any text in the red box that appears.



For further information on how to annotate proofs, click on the **Help** menu to reveal a list of further options:





**Preisliste für Sonderdrucke (2014)**

Die unten genannten Preise gelten nur für Sonderdrucke, die im Lauf des Jahres 2014 und vor dem Druck des Heftes bestellt werden. Mindestbestellmenge: 50 Exemplare. Lieferung: ca. 3 Wochen nach Erscheinen des Werkes, aus dem der Sonderdruck stammt.

Die Anfertigung von Sonderdrucken nach dem Druck einer Zeitschrift ist wesentlich teurer, da sie im Verlag und in der Druckerei zusätzliche Arbeitsgänge erfordert. Bitte bestellen Sie Ihre Sonderdrucke daher rechtzeitig und in ausreichender Stückzahl.

Bei Abnahme von mehr als 500 Exemplaren gelten Sonderpreise, die wir auf Anfrage mitteilen werden.

**Ganze Hefte können von Autoren zu einem ermäßigten Preis bestellt werden.**

Verpackung und Versandkosten (außer Luftpost) sind in den unten genannten Preisen enthalten. Die Mehrwertsteuer wird extra berechnet.

Sonderdrucke/Hefte werden per Schiff nach Übersee verschickt. Sollten Sie an einer schnelleren Beförderung per Luftpost interessiert sein, vermerken Sie dies auf beiliegendem Bestellformular. Bitte beachten Sie aber, dass Ihnen für diesen Versandweg die Kosten in Rechnung gestellt werden.

**Price List for Reprints (2014)**

The prices listed below are valid only for orders received in the course of 2014 and before the proofs pass for press. Minimum order is 50 copies. Delivery time will be approximately 3 weeks after the date of publication.

The production of reprints after a journal has been published is considerably more expensive since it requires extra operations on the publisher's and printer's side. Therefore, authors are requested to order reprints early and in sufficient numbers.

If more than 500 copies are ordered, special prices are available upon request.

**Single issues are available to authors at a reduced price.**

The prices include mailing and handling charges (with the exception of the additional costs incurred for air mail delivery).

All Wiley-VCH prices are exclusive of VAT.

Reprints and issues are shipped by surface. If you are interested in receiving them by airmail please indicate this on the accompanying order form; however, please be aware that the costs incurred are considerably higher (surcharge Euro 25.00).

Sonderdrucke/Reprints Umfang (Seiten) Size (pages)	Preis bei Abnahme von/Price for orders of (in Euro)					
	50 Expl./copies	100 Expl./copies	150 Expl./copies	200 Expl./copies	300 Expl./copies	500 Expl./copies
1– 4	345,—	395,—	425,—	445,—	548,—	752,—
5– 8	490,—	573,—	608,—	636,—	784,—	1077,—
9–12	640,—	739,—	786,—	824,—	1016,—	1396,—
13–16	780,—	900,—	958,—	1004,—	1237,—	1701,—
17–20	930,—	1070,—	1138,—	1196,—	1489,—	2022,—
je weitere 4 Seiten/ for every additional 4 pages	147,—	169,—	175,—	188,—	231,—	315,—
<b>Hefte/Issues</b>	1 Exemplar/1 copy: 25,— Euro		PDF (high resolution) 330,— Euro			

Postage and handling charges included. All Wiley-VCH prices are exclusive of VAT.

Prices are subject to change.