Check for updates

OPEN ACCESS

EDITED BY Dasiel Obregon, University of Guelph, Canada

REVIEWED BY

Harkirat Singh, Guru Angad Dev Veterinary and Animal Sciences University, India Sara Moutailler, Agence Nationale de Sécurité Sanitaire de l'Alimentation, de l'Environnement et du Travail (ANSES), France Márcia Cristina de Sena Oliveira, Embrapa Pecuária Sudeste, Brazil

*CORRESPONDENCE Juan Mosqueda joel.mosqueda@uag.mx

Joei.mosqueda@uaq.mx

SPECIALTY SECTION

This article was submitted to Veterinary Infectious Diseases, a section of the journal Frontiers in Veterinary Science

RECEIVED 28 September 2022 ACCEPTED 28 October 2022 PUBLISHED 11 November 2022

CITATION

Lizarazo-Zuluaga AP, Carvajal-Gamez BI, Wilkowsky S, Cravero S, Trangoni M and Mosqueda J (2022) Development and standardization of a Loop-mediated isothermal amplification (LAMP) test for the detection of *Babesia bigemina*. *Front. Vet. Sci.* 9:1056355. doi: 10.3389/fvets.2022.1056355

COPYRIGHT

© 2022 Lizarazo-Zuluaga, Carvajal-Gamez, Wilkowsky, Cravero, Trangoni and Mosqueda. This is an open-access article distributed under the terms of the Creative Commons Attribution License (CC BY). The use, distribution or reproduction in other forums is permitted, provided the original author(s) and the copyright owner(s) are credited and that the original publication in this journal is cited, in accordance with accepted academic practice. No use, distribution or reproduction is permitted which does not comply with these terms.

Development and standardization of a Loop-mediated isothermal amplification (LAMP) test for the detection of *Babesia bigemina*

Andrea P. Lizarazo-Zuluaga^{1,2}, Bertha I. Carvajal-Gamez^{1,3}, Silvina Wilkowsky⁴, Silvio Cravero⁴, Marcos Trangoni⁴ and Juan Mosqueda^{1,3*}

¹Immunology and Vaccines Laboratory, C. A. Facultad de Ciencias Naturales, Universidad Autónoma de Querétaro, Queretaro, Mexico, ²Maestria en Salud y Producción Animal Sustentable. Facultad de Ciencias Naturales, Universidad Autónoma de Querétaro, Querétaro, Mexico, ³C.A. Salud Animal y Microbiologia Ambiental. Facultad de Ciencias Naturales, Universidad Autonoma de Queretaro, Queretaro, Mexico, ⁴Instituto de Agrobiotecnología y Biología Molecular (IABIMO). Instituto Nacional de Tecnología Agropecuaria (INTA), Consejo Nacional de Investigaciones Científicas y Tecnológicas (CONICET), Hurlingham, Buenos Aires, Argentina

Bovine babesiosis is a tick-borne disease caused by protozoan parasites of the genus Babesia. Babesia bigemina is one of the most prevalent and economically important parasite species that infects cattle because of its impact on the meat and milk production industry. Effective disease control strategies should include detection of reservoir animals and early and specific pathogen detection using rapid, economical, sensitive, and specific detection techniques. The loop-mediated isothermal amplification technique (LAMP) is a one-step molecular reaction that amplifies DNA sequences with high sensitivity and specificity under isothermal conditions and requires no special equipment. The results can be observed by the naked eye as color changes. The aim of this work was to develop and standardize the LAMP technique for B. bigemina detection and its visualization using hydroxynaphtol blue. For this situation, primers were designed from the conserved sequences of the *B. bigemina ama-1* gene. The results showed that at 63 °C in 1 h and under standardized conditions, this technique could amplify B. bigemina DNA as indicated by the characteristic colorimetric change. Sensitivity evaluation indicated that DNA was amplified at a 0.00000001% parasitemia, and it was demonstrated that this technique specifically amplified the DNA of B. bigemina. Additionally, this technique could amplify DNA from 10 strains of *B. bigemina* from three different countries. It is concluded that the LAMP technique as modified in our case could specifically amplify *B. bigemina* DNA and shows high sensitivity, does not cross-react with related organisms, and the product is observed by 60 min of reaction time based on color changes. This report is the first LAMP report that uses sequences that are conserved between strains of the *ama-1* gene, demonstrates the results by color changes using hydroxynaphtol blue. We propose LAMP as a rapid and economical alternative method for the molecular detection of *B. bigemina*.

KEYWORDS

Babesia bigemina, molecular detection, LAMP, hydroxynaphtol blue, diagnostics

Introduction

Bovine babesiosis is one of the most economically important infectious diseases in the world. This disease is caused by intraerythrocytic parasitic protozoa of the genus Babesia (1, 2) and is transmitted by ticks of the genus Rhipicephalus (3, 4). Of the Babesia species that infect cattle, B. bovis and B. bigemina are the most prevalent and economically important, not only due to associated production losses and treatment costs but also because of trade restrictions and costs related to vector control (5). Therefore, effective disease control strategies must include detection of reservoir animals and early and specific diagnosis of the disease through rapid, sensitive, and specific detection techniques to establish treatment and/or control measures (2, 6-8). Most of the currently used techniques for the identification of B. bigemina have limitations: (1) the accuracy of the result is based on the training and expertise of the laboratory technician, (2) some methods need to be done during the acute phase of the disease to detect the parasite, and (4) several widely used techniques have low sensitivity and specificity (9, 10).

Molecular detection methods are characterized by high sensitivity and specificity for the detection of Babesia DNA in bovine blood; however, they are expensive techniques that require specialized equipment in addition to procedures that require time for the interpretation of a result. The loop-mediated isothermal amplification technique (LAMP), developed by Notomi et al. (11), is a relatively simple method that does not require specialized equipment because the amplification is isothermal and requires only a single step, and in 45 min, it amplifies the target DNA up to 109 copies with high specificity, efficiency, and speed. The LAMP method uses the Bst DNA polymerase (an enzyme that has chain displacement activity) and six to eight oligonucleotides called external oligonucleotides, which provide the technique with specificity and sensitivity; additionally, it uses LOOP oligonucleotides (LP and LB), a step that optimizes amplification time (12). Detection of the amplified product can be confirmed by visualization on agarose gels, by fluorescent dyes, or with chemical agents, such as hydroxynaphthol blue (HNB), which produces a colorimetric change (13, 14). HNB is a metal indicator that changes color from violet in a negative reaction to blue in a positive reaction resulting from changes in magnesium's ion (Mg) concentration.

Several advantages of using HNB should be mentioned: (1) it does not require opening the reaction tube; therefore, it is not necessary to analyze it by electrophoresis and (2) the result can be differentiated with the naked eye (13, 14). An important element for the success of the LAMP technique is the selection of the target sequence and the design of specific oligonucleotides (14). Previous studies have reported the development of a LAMP technique for the detection of B. bigemina and B. bovis. In these studies, they used sequences, such as the gene that encodes the rhoptry-associated protein 1 (rap-1) and an intergenic internal spacer sequence (ITS), which are characteristic of variable regions of ITS sequences. These techniques were evaluated in bovine samples located in different locations in Ghana and Zambia in addition to China (10, 15). Additionally, a LAMPlateral flow dipstick (LAMP-LFD) technique was developed to detect the gene that encodes cytochrome b in *Babesia* spp. (16). The gene that encodes the merozoite apical antigen 1 (ama-1) is a gene whose product is expressed on the surface of merozoites and is responsible for interacting with the rhoptry neck protein (RON2) protein and forming a tight junction between Babesia and an erythrocyte. Studies have shown that the gene sequence presents variable or polymorphic regions, observed mainly in domains I and III, and conserved regions in domain II. Additionally, these studies have shown its high level of conservation and low level of polymorphisms (17-19). The ama-1 gene is the proposed target in this work due to its ideal properties for use in molecular detection techniques, such as containing specific sequences of B. bigemina and its conservation among all strains sequenced to date (17).

The aim of this study was to develop and standardize a LAMP-based technique for molecular detection of *B. bigemina* using *ama-1* as a target in addition to its visualization using HNB.

Materials and methods

Bioethical declarations

The protocol for handling and bleeding the cattle and for tick collection was approved by the Bioethics Committee of the College of Natural Sciences, Autonomous University of Queretaro, Mexico (approval numbers; FCN/2011-0221 and 71FCN2016).

Biological samples

To obtain the DNA of B. bigemina, a male Holstein calf of approximately 1 year of age from a tick-free area and negative for the presence of antibodies against B. bovis, B babesia, B. bigemina, and A. marginale was splenectomized and experimentally infected with a stabilate containing B. bigeminainfected erythrocytes, Chiapas strain, which had been stored in liquid nitrogen. When the calf reached 2% parasitaemia, blood was obtained in a sterile flask with glass beads. Blood was washed three times via centrifugation with VYM solution to remove white blood cells, and the infected erythrocytes were stored at -20 °C until use. Additionally, Rhipicephalus microplus ticks infected with Babesia spp. were collected from the states of Chiapas, Colima, Jalisco, Nayarit, and Queretaro in Mexico. The eggs of the infected ticks were incubated at 28 °C until the larvae hatched. The larvae were placed on a calf, and 6 days later the engorged larvae were removed and incubated in the laboratory at 28 °C until they molted into nymphs. The nymphs of each isolate were placed in a splenectomized calf to transmit B. bigemina. When intraerythrocytic merozoites were observed, blood was collected and processed in the same way as described above for the Chiapas strain. Each calf received anti-babesia treatment after blood collection. On the other hand, genomic DNA samples of the S1A, S2P, and S3P strains from Argentina, which had been stored frozen in the Hemoparasites Laboratory of the Institute of Agrobiotechnology and Molecular Biology (INTA-CONICET), in Argentina, were used. Finally, an isolate of B. bigemina obtained from a steer in Colombia that had clinical babesiosis was also included in this study.

Genomic DNA

Total DNA extraction was performed using the Illustra Blood Genomic Prep mini spin kitTM (GE Health Care Life Sciences, USA) following the manufacturer's instructions, using 200 μ l of blood. Genomic DNA was quantified in a NanoDrop 2000TM (Thermo Fischer Scientific, Waltham, MA, USA) equipment. DNA with purity between 1.7 and 1.8, which is the value of the ratio of the absorbance reading measured at 260 and 280 nm, was selected and stored in small aliquots at -20 °C until use.

Bioinformatic analysis of the *ama-1* gene and oligonucleotide design

To design *B. bigemina*-specific oligonucleotides, the single copy *ama-1* gene was selected because it has a low number of polymorphisms and a high level of conservation among the different strains of *B. bigemina* (17). Subsequently, a BLAST (NCBI) (https://blast.ncbi.nlm.nih.gov/Blast.cgi) was performed

using the sequence of the ama-1 gene from the Chiapas and Australian strains, which are available on the Wellcome Trust Sanger Institute website (https://www.sanger.ac.uk/resources/ downloads/protozoa/babesia-bigemina.html). A search of sequences reported to date was performed in the GenBank (20) available online at www.ncbi.nlm.nih.gov/genbank/. Thirtyone sequences were used on which multiple alignment was performed with the Clustal Omega program of the European Institute of Bioinformatics (https://www.ebi.ac.uk/Tools/ msa/clustalo/). From the multiple alignment, the consensus sequence, which includes the conserved regions of the ama-1 gene, was determined. This sequence was analyzed using the PrimerExplorer program, version 5.0 (21). Subsequently, each set of resulting oligonucleotides was evaluated using the OligoAnalyzer program, version 3.1 from Integrated DNA Technologies (IDT), and it was verified that the set met the specifications reported by (12) and the specificity of each oligonucleotide (22). Selected oligonucleotides were synthesized commercially (Genbiotech, Buenos Aires, Argentina), high-performance liquid chromatography (HPLC) grade.

LAMP reaction standardization and optimization

To perform the standardization of the LAMP method, DNA from B. bigemina, Chiapas strain, was used. The reaction mixture consisted of 1.6 µM inner primers (FIP and BIP), 0.2 µM outer primers (F3 and B3) and 0.8 µM loop primers (FLP and RLP), a final concentration of 1.4 mM of deoxyribonucleotides triphosphate (dNTPs) from Invitrogen, Life Technologies (MA, USA), 0.8 M betaine (Sigma-Aldrich, Burlington, MA, United States), 1x buffer (New England Biolabs, Massachusetts, USA), 6 mM MgSO4 (New England Biolabs), 8U of Bst polymerase (New England Biolabs. MA, USA), 120 µM of HNB (Sigma-Aldrich, Burlington, MA, USA), and 1 µL containing 100 ng of DNA and nuclease-free water to complete 25 µL final volume. The LAMP reaction was evaluated at temperatures between 60 and 65 °C, at different incubation times (15, 30, 45, and 60 min) with a final incubation temperature of 83 °C for 3 min. Additionally, different concentrations of oligonucleotides (ratio of internal to external primers at an external:internal ratio of 1:2, 1:3, 1:4, 1:5, and 1:6) were evaluated.

Specificity and sensitivity determination

Once the optimization was done, diagnostic specificity and sensitivity were determined with a reaction mixture consisting of $1.6 \,\mu$ M inner primers (FIP and BIP), $0.2 \,\mu$ M outer primers (F3 and B3) and $0.8 \,\mu$ M loop primers (FLP and RLP) at a final concentration of the same reagents described above and the selected target DNA. The diagnostic sensitivity of the technique was determined by performing serial dilutions of the DNA extracted from a blood sample with 2% parasitemia using dilutions in a range of 10^{-2} - 10^{-9} . Additionally, the diagnostic specificity of the LAMP technique was performed using DNA from *B. bovis, Anaplasma marginale, A. phagocytophilum, A. centrale, Trypanosoma theileri, Bos taurus, Homo sapiens, Rhipicephalus microplus* and *Neospora caninum.* The results were determined visually according to the procedures reported by Goto et al. (13) using agarose gels stained with GelRed (Sigma-Aldrich, Burlington, MA, United States) and analyzed using a gel documentation equipment (Bio-Rad. Hercules, California, United States).

Comparison of LAMP with the gold standard detection method

To compare the sensitivity of the LAMP developed in this study with the nested PCR (nPCR), which was considered the gold standard for the molecular diagnosis of *B. bigemina*, the primers, procedures, and programs reported by the World Organization for Animal Health, (23) and by Figueroa et al. (24) were used as references. The analysis of the results was done visually (13) and on an agarose gel stained with GelRed using a gel documentation equipment (Bio-Rad. Hercules, California, United States).

Evaluation of the LAMP test for the detection of *B. bigemina* strains from different geographical areas

Once the technique was optimized, it was used to detect DNA of *B. bigemina* strains from different regions of the world. The experiment was performed with DNA from strains from Mexico, including Chiapas, Colima, Jalisco, Nayarit and Queretaro, strains from Argentina, which included S1A, S2P and S3P, and DNA of an isolate from Colombia.

Results

Primer design

The oligonucleotides proposed in this work were designed based on a conserved sequence of the *ama-1* gene. The multiple alignment used 31 sequences from different strains of *B. bigemina* and showed a high degree of conservation of 100% between the sequences from different areas of the world (Figure 1). The selection of the oligonucleotides was shown to be species-specific using BLAST-based amplification of a 204 bp fragment of the *ama-1* gene of *B. bigemina*. Additionally, the proposed oligonucleotides fulfilled the characteristics required for the LAMP test. The sequences of the oligonucleotides are shown in Table 1.

Standardization and optimization of the LAMP reaction

The standardization of the technique using total DNA from blood samples with a 2% parasitemia of B. bigemina showed amplification at the six experimental temperatures. This result shows that the LAMP test can be performed in this temperature range in a water bath. Additionally, visualization showed colorimetric change at all experimental temperatures, which was represented as a change in color from violet (negative reaction) to sky blue (positive reaction) without differences between each one. The result agreed with the electrophoretic image (Figures 2A,B). From this point, the reactions were incubated at 63 °C, as the optimal temperature. Regarding the minimum time necessary for DNA amplification using the LAMP technique, the result shows that after 15 min of reaction, no amplification had occurred. However, at 30 min, the characteristic amplicons of the LAMP reaction were observed (Figure 3A). Additionally, the evaluation of the colorimetric change was observed from 30 min as a color change from violet to blue (Figure 3B).

Evaluation of the specificity and sensitivity of the LAMP test

The specificity of the test showed that the oligonucleotides are specific for B. bigemina since no amplification or colorimetric changes were observed in the reactions using the DNA of the evaluated microorganisms (Figures 4A,B). This result agrees with the results obtained by bioinformatic analyses. On the other hand, the sensitivity of the test when evaluating the dilutions using DNA extracted from blood containing 2% parasitaemia showed a detection limit of 1×10^{-8} in the LAMP test, indicating that it can detect 0.00000001% parasitemia in whole blood samples (Figure 5A). Importantly, a gradual color change from light blue to dark blue was observed according to the dilution evaluated visually (Figure 5B). Additionally, the 1×10^{-9} dilution showed a violet coloration, indicating a negative amplification reaction. The results obtained by electrophoresis and those obtained visually agree in both experiments (Figures 5A,B), which suggests that both visualization methods have the capacity to detect the amplified product in the same dilution. Therefore, the results show that both techniques possess the same sensitivity for amplification detection. In addition, when comparing the LAMP and the nPCR techniques, it was shown that both techniques

Australian	GASCITCSSCASSIGCTCTTACAAACSCASSCASCATCCACSSIGAATTTGTOSCCASTAT CASCOSCOSACCTCCSTAASTSSSSTACAASSSAAACSCCSTTACTAACTSSCCSASTACTSCCACSAATACSCCSASTACTSCCASSAACSCCSACACSCCSACACSCCSACACSCCSACACSCCSACCTCCSTAASTSSSSTACSCCSACACSCCSACCCCSACSCCSACCTCCSTAASTSSSSTACSCCSACACSCCSACCCCSACSCCSACCSCCSACSCCSACCSCCS
SP3	GAGCTTC09CA9GT0CTCTTACAAAC9CA9GCA9CATCCAC9GTAATTTGTC9CCAGTAT CA9CCDCC6ACCTCC6TAAST6995TTACAA956AAAC9CC9TTACTAACT6C9CC6AATAC9CCA9CA
172	GASCITCSSCASSIGCTCTTACAAACSCASSCAACATCCACSSIAATTTSTCSCCASIAT CASCSCCSACCITCCSIAASTSSSSTACAASSCSCCSTIACTAACTSCCCCSAATACSCCSACA
IV1	GAGCTTC09CA9GT0CTCTTACAAAC9CA9GCAACATCCAC9GTAATTTGTC9CCAGTAT CA9CC9CC6ACTCCGTAAGT099STTACAA99GAAAC9CCGTTACTAACT9C9CCGAAT AC9CCA9CA
ITASCT	GASCTTC05CASGTGCTCTTACAAAC9CAACATCCAC05TAATT0FC05CCASTAT CA9CC9CC05AACT0505TAA5T99095TTACAA6953AAC9CC9TTACTAACT9C9CC5GAATAC9CCA9CA
ITA7PA	GAGCTTC09CA9GT0CTCTTACAAAC9CCA9GCAACATCCAC9GTAATTTGTC9CCAGTAT CA9CCSCCSACCTCCGTAAGT909STTACAA9SGAAAC9CCGTTACTAACT9C5CCCGAAT AC9CCA9CA
ITA6PA	GAGCTTCGSCAGGTGCTCTTACAAACGCAGGCAACATCCACGGTAATTTGTCGCCAGTAT CAGCCGCCGAACTCCGTAAGTGGGGTTACAAGGGGAAACGCCGGTAACTGCGCCGCGAAT ACGCCAGCA
ITA58G	GAGCTTCGSCAGGTGCTCTTACAAACGCAGGCAACATCCACGGTAACTTGTCGCCAGTAT CAGCCGCCGAACTCCGTAAGTGGGGTTACAAGSGAAACGCCGTTACTAACTGCGCCGCGAAT ACGCCAGCA
ITA4PG	GAGCTTCGSCAGGTGCTCTTACAAACGCAGGCAACATCCACGGTAATTTGTCGCCAGTAT CAGCCGCCGAACTCCGTAAGTGGGGTTACAAGGGGAAACGCCGGTTACTAACTGCGCCGCGAAT ACGCCGCAGTAT
Turkey	GAGCTTCGSCAGGTGCTCTTACAAACGCAGGCAACATCCACGGTAATTTGTCGCCAGTAT CAGCCSCCGAACTCCGTAAGTGGGSTTACAAGSGAAACGCCGTTACTAACTGCGCCGCGAAT ACGCCAGCA
Rio	GASCTTCOSCASSTGCTCTTACAAACSCASGCAACATCCACOSTAATTGTCOSCCASTAT CASCCSCCOSCACCTCCSTAASTGSSSTACAASSSAAACSCCSTTACTAACTGCSCCSAAT ACSCCASCA
Mexico	GASCITCOSCASSIGCTCITACAAACSCASCAACATCCACOSTAATITGTOSCCASTAT CASCOSCOSACCTCCSTAASTSSSSTACAASSSAAACSCCSTTACTAACTSCSCCSAATACSCCASCA
MIA	GASCTTCSSCASSTGCTCTTACAAACSCASCAACATCCACGSTAATTGTCSCCASTAT CASCCSCCSACTCCGTAASTG995TTACAAS953AACSCCGTTACTAACT9C5CCSAAT AC5CCASCA
M1P	GAGCTTC09CA9GT9CTCTTACAAAC9CCA9GCAACATCCAC9GTAATTTGTC9CCAGTAT CA9CCSCCSACCTCC5TAASTG995TTACAA9563AAAC9CCSTTACTAACT9C5CCC5AAT AC9CCA9CA
ITA3	GASCITCSSCASSIGCTCTTACAAACSCASSCAACATCCACSSIGAATTIGTCSCCASTAT CASCSSCCSACCTCCSTAATSGSSGTACAASSSAAACSCCSTTACTAACTSSCCSSATACSCCSACA
ITA2	GASCITCSSCASSTSCICTTACAAACSCASSCAACATCCACSSGTAATITGTCSCCASSTAT CASCCSCCSACCTCCSTAASTOSSSTTACAASSSAAACSCCSTTACTAACTSCCCCSGAAT ACSCCASCAAT ACSCCASCAAT
ITA1	GAGCTTCSGCAGGTGCTCTTACAAACSCAGGCAACATCCACSGTAATTTGTCGCCAGTAT CAGCCSCCSACCTCCSTAASTG9995TTACAA999AAAC9CCSTTACTAACT9CSCCSAAT AC90CA9CA
Bb1MF-C4	GAGCTTCGGCAGGTGCTCTTACAAACGCAGGCAACATCCACGGTAATTTGTCGCCCAGTAT CAGCCGCCGGACCTCCGTAAGTG95GTTACAAG9GAAACGCCGTTACTAACTGC5CCGAAT AC9CCAGCA
GP-C19	GASCITCSSCASSISCECTIACAAACSCASSCAACATCCACSSTAATTTSTCSCCASTAT CASCCSCCSACTCCSTAAST909STTACAASSSAAACSCCSTTACTAACT90SCCSAAT A090CASCA
BbiGP-C7	GASCITCOSCASSIGCTCTTACAAACSCASSCAACATCCACSSTAATTIGTOSCCASTAT CASCOSCOSACCTCCSTAASISSSAAACSCOCSTTACTAACTSCSCCSAAT ACSCCASCA
GP-C1	6ASCITCSSCASSISCICTIACAAACSCASSCAACATCCACSSIAATITSTOSCCASIAT CASCOSCOSACITCCSIAASISSSAAACSCOSITACAASSCAACACSCOSTAATACSCOSCOSAATACSCOSCASAA
FS-189	GASCITCOSCASSIGCTCTTACAAACSCASSCAACATCCACSSIAATTTGTCSCCASTAT CASCCSCCSACCTCCSTAATSGSSITACAASSCSAACSCCCGTTACTAACTSCSCCSAAT ACSCCASCAJ
FS-146	GASCITCOSCASSIGCTCTTACAAACSCASSCAACATCCACSSIAATTTSTCSCCASTAT CASCOCCCCSACCTCCSTAASIGSSATACSCCCCTTACTAACTGCSCCSAAT ACSCCASCA
WC-10284	GAGCTTC09CA9GT0CTCTTACAAAC9CCA9GCAACATCCACGGTAATTTGTC9CCAGTAT CA9CC9CC9ACTCCGTAAGT099STTACAA9SGAAAC9CCCGTTACTAACT9C9CCCGAAT AC9CCA9CA
WC-851	GASCITCOSCASOTOCICITACAAACOCAOOCAACATCCACOOTAATITOTCOCCAGTAT CASCOCCCGACCICCGTAAGTOGOGITACAAGOGAAACOCCGTTACTAACTOCCCCGAATACCOCCGACACACCACCGACACACCACCGACACACCCCCGACACACCCCGACACACCCCGACACACCCCCGACACACCCCCGACACACCCCCGACCACC
MP-C12	GASCITCOSCASGIGCICCITACAAACOCAOGCAACATCCACOGTAATITGICOCCAGIAT CASCOCCCGACCITCCGIAAGIGOSOTIACAAOSOCGATACGCCGITACIAACIOCCCCGAAI ACOCCAOCA
MP-C11	GASCITCOSCASSIGCTCITACAAACSCASSCAACATCCACSSIAATITGTCSCCASIAT CASCOCCCGACTCCGTAASIGSSITACAASSGAAACSCCGTTACTAACTGCSCCGAATACSCCASCA
MP-C1	GASCITCSSCASSIGCTCITACAAACSCASSCAACATCCACSSIAATITGTCSCCASIAT CASCOSCCSACCTCCSIAASISSSSTACAASSCSCCSITACIAACISCCSCCSAAIACSCCSICAAAACSCCSICAAAACSCCSICAAAACSCCSICAAAAACSCCSICAAAAACSCCSICAAAAACSCCSICAAAAACSCCSICAAAAACSCCSICAAAAACSCCSICAAAAACSCCSICAAAAACSCCSICAAAAACSCCSICAAAAACSCCSICAAAAAACSCCSICAAAAAACSCCSICAAAAAACSCCSICAAAAACSCCSICAAAAAACSCCSICAAAAAACSCCSICAAAAAACSCCSICAAAAAACSCCSICAAAAAACSCCSICAAAAAAACSCCSICAAAAAAACSCCSICAAAAAACSCCSICAAAAAACSCCSICAAAAACSCCSICAAAACSCCSICAAAAAACSCCSICAAAAAACSCCSICAAAAAACSCCSICAAAAACSCCSICAAAAAAAA
DIVA	GASCITCSSCASSISCITCITACAAACSCASSCAACATCCACOSTAATITGTCSCCASTAT CASCOSCCSACCITCCSTAASISSSSTACAASSCAAAACSCCSTTACTAACIOCSCCSAATACSCASCA
KZN-C24	GASCITCSSCASSTSCTCTTACAAACSCASSCAACATCCACSSTAATTIGTCSCCASTAT CASCCSCCSACCTCCSTAAST9999TTACAA993AAACSCCSTTACTAACT9CSCCSSATACSCASCA
Corrientes	GASCITCSSCASSIGCTCTTACAAACSCASSCAACATCCACSSITAATITISTCSCCASSIAT CASCCSCCSACCTCCSITAASIGSSSITACAASSCSAAACSCCSITIGCIAACTSCCSAATACSCCSSAATACSCCSSCA
Chispas	ATCGTGCCCGGTTCCGACCAGBGACCAAGTACAGGTACCCCTTCGTCT ATGATGGAAGGAGGAGATGTGTTACATCCTGTACTCGCCGATGCAGTACCAGGGAACCAGGGAACCAGGGAACCAGGGAACCAGGGAGACCAGGGAGACCAGGAGG
Australian	A TO STOCCOSTICCONCEANSONCEAN STACKOGIA CONCISCION OF A GALGGAA SGAGGAGA STOFFIACA CONSTACT STACKOGGA CASSA
SP3	A TOSTOCCOSSI TO GALCASASSACASSACASSACCOST COST A CAST GALASSASSASATO STILLATO TO TACTOSCOST SCASSACASSAA CASSIACTO TO ACCASSAC
112	A TO STOCCOSTICCONCLAMAGINE CAMPACAMOTACCONTONE TA CAN SO A CAN STACT STACAMONA SO A S
IV1	AICOTGCCC0GTCCC3AACAAGCACAAGTACAGGTACCCCTTCGTCT ACGATGGAAAGGAGGAGATGTGTIACATCCTGTACTCGCCGATGCAGTACCAAGCGAAGAACCAGGGAACCAGGGAACCAGGGACGAGGAC
ITASCT	ATCSTGCCCGGTTCCGACCAGAGGACCAAGTACAGGTACCCCTTCGTCT ACGATGGAAAGGAGAGATGTGTTACATCCTGTACTCGCCGATGCAGTACCAGGAACCAGGTACTGTGACGAGGAG
ITA7PA	ATCSTGCCCGGTTCCSACCAGAGGACCAAGTACAGGTACCCCTTCGTCT ACGATGGAAGGGGGGGGGTGCGTTACATCCTGTACCGGCCGATGCAGTACAACCAGGGAACCAGGTACTGTGACGAGGACGAGGAC
ITA6PA	
	ATCSTSCCC8STTCC8ACCASSACCASSACCASSTACCCCTTCSTCT AUXAL 90AAA99A99A94191911AUA10UU191AU10UU419UA91AUA8UUA994AACCASSTACTSTGACSASSAC
ITASRG	
	ATCSTGCCCGGTTCCSACCAGAGGACCAAGTACAGGTACCCCTTCGTCT ACGATGGAAAGGAGGAGATGTGTTACATCCTGTACTCGCCGATGCAGTACAACCAGGGAACCAGGTACTGTGACGAGGAGAC
ITASRG	ATCSTGCCC93TTCC3ACCA393SACCAASTACA39TACC3CTTCSTCT ACSATGSAAA93A93A93A93A93A9TGTTTACATCCTGTACTC9CC3AT9CA3TACAA9CA39SAACCA3
ITA58G ITA4PG	ATOSTOCCOSTTCOSACCASSSACCASTACASSTACCOSTCOSTCT ACOATOSAAASASAGASTGTTACATCCTGTACTOSCOSTCOSTCASTACAACASOSAA CCASSTACTSGTAACASOSAC ATOSTOCCOSSTTCOSACCASSSACCASTACASSTACASSTACCOCTTOSTCT ACOATOSAAASASASAGASTGTTACATCCTGTACTOSCOSTGCASTACAACASOSAA ATOSTOCCOSSTTCOSACCAASSACCAASTACASSTACCOCTTOSTCT ACOATOSAAASASASAGASTGTTACATCCTGTACTOSCOSTGCASTACAACASOSAA CCASSTACTGTGACGASSACC
ITASRG ITA4PG Turkey Rio	ATCSTGCCCGSTTCCSACCASSGACCAASTACAGGTACCCCTTCGTCT ACGATGGAAAGGAGGAGATGTGTIACATCCTGTACTCGCCGATGCAGTACAACCAGGGAACCAGGTACTGTGACGGGGAG ATCSTGCCCGGTTCCGACCAASGAACCAASTACAGGTACCCCTTCGTCT ACGATGGAAAGGAGGAGATGTGTIACATCCTGTACTCGCCGATGCAGTACCACGGGAACCAGGGAACCAGGTACCGCGGTCCGACCAGGGAGACCAAGTACAGGTACCCTGGTCT ACGATGGAAGGGGGGAGATGTGTIACATCCTGTACTCGCCGATGCAGTACAACGAGGAACCAGGGAACCAGGTACCGTGGACGAGGACGAGGACGAGGGGGGGG
ITASRG ITA4PG Turkey Rio Mexico MlA	ATOSTOCCOSTICOLACCABASAACAASTACASOTACCOCTTOSTCT ACGATOSAAASAASAATSTOTTIACATCOTOTACTOSCOATGCASTACAACCASOSAA CCASOTACTOSTAACGASAAC ATOSTOCCOSTICOLACCABASGACCAASTACASOTACCOCTTOSTCT ACGATOSAAASAAGAGAGATGTSTTACATCOTOSCOATGCASTACCAACGASGAACCASOTACTOSTGACGASGAC ATOSTOCCOSTICOLACCAASSAACCAASTACASOTACCOCTTOSTCT ACGATOSAAASGAGAGATGTSTTACATCOTOTACTOSCOATGCASTACCAACGASGAACCASOTACTOSTGACGASGAC ATOSTOCCOSTICOLACCAASSAACCAASTACASOTACCOCTTOSTCT ACGATOSAAASGAGAGATGTSTTACATCOTGTACTOSCOATGCASTACCAACGASGAACCASOTACTGTGACGASGAC ATOSTOCCOSTICOLACCAASSGACCAASTACASOTACCOCTTOSTCT ACGATOSAAASGAGAGAGATGTSTTACATCOTGTACTOSCOATGCASTACAACCASOSAACCASOTACTGTGACGASGAC ATOSTOCCOSTICOLACCAASGACCAASTACASOTACCOCTTOSTCT ACGATOSAAASGAGAGAGATGTSTTACATCOTGTACTOSCOATGCASTACCASGAACCASOSAACCASOTACTGTGACGASGAC ATOSTOCCOSTICOLACCAASGACCAASTACASOTACCOCTTOSTCT ACGATOSAAGGAGAGAGAGATGTSTTACATCOTGTACTOSCOATGCASTACCASGGAACCASOTACTGTGACGASGAC ATOSTOCCOSTICOLACCAASGACCAASTACASOTACCOCTTOSTCT ACGATOSAAGGAGAGAGATGTSTACATCOTGTACTOSCOATGCASTACCASGGAACCASGAACCACGASGACCAGCACGASTACTGTGACGASGACGAGAGAGAGAGAGAGAGAGAGAGAGAGAGAG
ITASRG ITA4PG Turkey Rio Mexico MIA MIP	ATCOTOCCOGNTCOSACCARGORACCARGTACABOTACCCCTTODICT ACORTOGRAMAGRAGORATOTGTIACATCCTOTACTCOCCORTOCAGTACAACCAGOGRA CCARGTACTGTGACGAGAAC ATCOTOCCOGNTCOSACCARGORACCARGTACABOTACCCCTTOGICT ACORTOGRAMAGRAGARGTOTIACATCCTOTACTOSCCORTOCAGTACAACCAGOGRA CCARGTACTGTGACGAGAAC ATCOTOCCOGNTCOSACCARGORACCARGTACABOTACCCCTTOGICT ACORTOGRAMAGRAGARGTOTIACATCCTOTACTOSCCORTOCAGTACAACCAGOGRA CCARGTACTGTGACGAGAAC ATCOTOCCOGNTCOSACCARGORACCARGTACABOTACCCCTTOGICT ACORTOGRAMAGRAGARGTOTIACATCCTOTACTOSCCORTOCAGTACAACCAGOGRA CCARGTACTGTGACGAGAAC ATCOTOCCOGNTCOSACCARGORACCARGTACABOTACCCCTTOGICT ACORTOGRAMAGRAGARGARGTOTIACATCCTGTACTOSCCORTOCAGTACAACCAGOGRA CCARGTACTGTGACGAGAAC ATCOTOCCCGNTTCOSACCARGORACCARGTACABOTACCCCTTOGICT ACORTOGRAMAGRAGARGARGARGTOTIACATCCTGTACTOSCCORTGCARGTACAACCABOGRACCAGGAACCABOGAACTGTGTACATCCTGTACACCACGAGAACCARGGARGTACTGTGAGAGAACAACGAGAACCAGOGRA CCARGTACTGTGACGAGAACCARGTACTGTGACGAGAACCARGTACTGTGACGAGAACCARGGAACTGTGTACATCCTGTACCACGTACCACGAGAACCARGGAACCARGGAACCAGGAACCAGGAACCARGTACTGTGACGAGAACGAGAACGACGAACGACGAACGACGACCAACGAGAACCARGGAACGACGACGACGACGACGACGACGACGACGACCARGTACCACGAGAACGACGACGACGACGACGACGACGACGACGA
ITASRG ITA4PG Turkey Rio Mexico MIA MIP ITA3	ATCSTGCCCGSTTCCSACCAGAGGACCAASTACAGGTACCCCTTCGTCT ACGATGGAAAGGAGGAGATGTGTTACATCCTGTACTCGCCGATGCAGTACAACCAGGGAACCAGGTACTGTGACGAGGAC ATCGTGCCCGGTTCCGACCAGAGGACCAASTACAGGTACCCCTCGTCT ACGATGGAAAGGAGGAGGAGGAGGAGTGTTACATCCTGTACTCGCCAGTACAACCAGGGAA ATCGTGCCCGGTTCCGACCAGAGGACCAASTACAGGTACCCCTTCGTCT ACGATGGAAAGGAGGAGGAGGAGGAGTGTTACATCCTGTACTCGCCGAGTACAGACGAGGACCAGGTACTGTGAGGAGGAC ATCGTGCCCGGTTCCGACCAGAGGACCAASTACAGGTACCCCCTTCGTCT ACGATGGAAAGGAGGAGGAGGAGGAGGAGGAGGAGGAGGAGG
ITASRG ITA4PG Turkey Rio Mexico MiA MiP ITA3 ITA2	ATOSTOCCOSTTCOLACCABASIACCASTACASSTACCOSTCOTT ACIATOSIAASIASIAATOTTIACATCITACATCITACATCASTACAACCASSBACCASSTACTSTAACGASSAC ATOSTOCCOSTTCOLACCABASIACCASTACASSTACCOSTCOTTCACIATOSIAASIASIASIASTOTTIACATCITACATCITOCCOATOCASTACAACCASSBAACCASSTACTSTGAGASSAC ATOSTOCCOSTTCOLACCASSIACCASTACASSTACCOCTTOSTCT ACIATOSIAASIASIASIASTOTTIACATCITGTACTOSCOATOCASTACAACCASSBAACCASSTACTSTGAGASSAC ATOSTOCCOSTTCOLACCASSIACCASTACASSTACCOCTTOSTCT ACIATOSIAASIASIASIASTOTTIACATCITGTACTOSCOATOCASTACAACCASSBAACCASSTACTSTGAGASSAC ATOSTOCCOSSTTCOLACCASSIACCASSTACASSTACCOCTTOSTCT ACIATOSIAASIASIASIASIASTOTTIACATCITGTACTOSCOATISCASTACAACCASSBAACCASSTACTSTGAGASSAC ATOSTOCCOSSTTCOLACCASSIASIACASSTACASSTACCOCTTOSTCT ACIATOSIAASIASIASIASTISTITACATCITGTACTOSCOATISCASTACAACCASSBAACCASSTACTSTGAGASSAC ATOSTOCCOSSTTCOLACCASSIASICAASSTACASSTACCOCTTOSTCT ACIATOSIAASIASIASIASTISTITACATCCTGTACTOSCOATISCASTACAACCASSBAACCASSBAACCASSIAS ATOSTOCCOSSTTCOLACCASSIASICAASSTACASSTACCOCTTOSTCT ACIATOSIAASIASIASTISTITACATCCTGTACTOSCOATISCASTACAACCASSBAACCASSACCASSTACTSTGAGASSAC ATOSTOCCOSSTTCOLACCASSBACCASSTACASSTACCOCTTOSTCT ACIATOSIAASIASIASIASTISTITACATCCTGTACTOSCOATISCASTACAACCASSBAACCASSBAACCASSTACTGTGAGASSAC ATOSTOCCOSSTTCOLACCASSBACCASSTACASSTACCOCTTOSTCT ACIATOSIAASIASIASIASTISTITACATCCTGTACCOCCUSTICCASTACCAACCASSBAACCASSTACTGTGAGASSAC ATOSTOCCOSSTTCOLACCASSBACCASSTACASSTACCOCTTOSTCT ACIATOSIAASIASIASIASTACTGTTACATCCTGTACCOCCUSTICCASTACCAACCASSBAACCASSBACCASSTACTGTGAGASSACCA
ITASRG ITASRG Turkey Rico MEXICO MIA MIP ITA3 ITA2 ITA1	ATCOTOCCOGNITCOSACCAASGAACAASTACABOTACCCCTTOGICT ACGATOGAAAGAGAGATOTGTIACATCCTGTACTCGCCGATOCAGTACAACCAGGGAACCAGGTACTGTGAGGAGAGAC ATCOTOCCOGNITCOSACCAASGAACCAASTACABOTACCCCTTOGICT ACGATOGAAAGAGAGATOTGTIACATCCTGTACTGCCCGATOCAGTACAACCAGGGAACCAGGTACTGTGAGGAGAGAC ATCOTOCCCGGTTCOSACCAASGAACCAASTACABOTACCCCTTOGICT ACGATOGAAAGAGAGATOTGTIACATCCTGTACTGCCCGATOCAGTACAACCAGGGAACCAGGTACTGTGAGAGAGAC ATCOTOCCCGGTTCOSACCAASGAACCAASTACABOTACCCCTTOGICT ACGATOGAAAGAGAGATGTGTIACATCCTGTACTGCCCGATOCAGTACAACCAGGGAACCAGGTACTGTGAGAGAGAC ATCOTOCCCGGTTCOSACCAASGAACCAASTACABOTACCCCTTOGICT ACGATOGAAAGAGAGATGTGTIACATCCTGTACTGCCCGATOCAGTACAACCAGGGAACCAAGGAACCAGGTACTGTGAGAGAGA
ITASRG ITA4PG Turkey Rio Mexico MiA MiP ITA3 ITA2	ATOSTOCCOSTITCOLACCABASIACCASTACASOTACCOCTTOSTCT ACIATOSIAACSASGASTITOTTACATCOTOTACTOSCOATOCATICASTACAACCASOSIAACCASOTACCOCTOSTCT ACIATOSIAACSASGASTACTOTTACATCOTOTACTOSCOATOCASTACAACCASOSIAACCASOTACTOTGASGASGAS ATOSTOCCOSTITCOLACCASSSACCASTACASOTACCOCTTOSTCT ACIATOSIAASGASGASTOSTITACATCOTOTACTOSCOATOCASTACAACCASOSIAACCASOTACTOTGASGASGA ATOSTOCCOSTITCOLACCASSSACCASTACASOTACCOCTTOSTCT ACIATOSIAASGASGASTOSTITACATCOTOTACTOSCOATOCASTACAACCASOSIAACCASOTACTOTGACGASGAC ATOSTOCCOSTITCOLACCASSSACCASTACASOTACCOCTTOSTCT ACIATOSIAASGASGASGASTOSTITACATCOTOTACTOCCOSTICOLASTACASCASSGACCASSTACTOTGACGASGAC ATOSTOCCOSTITCOLACCASSGACCASSTACASOTACCOCTTOSTCT ACIATOSIAASGASGASGASTOSTITACATCOTOTACTOCOCOSTICOLASTACASCASSGACCASSTACTOTGACGASGAC ATOSTOCCOSTITCOLACCASSGACCASSTACASOTACCOCTTOSTCT ACIATOSIAASGASGASGASTOSTITACATCOTGACTOSCOATOCASTACASCASSGAACCASSTACTOTGACGASGAC ATOSTOCCOSSTICOLACCASSGACCASSTACASOTACCOCTTOSTCT ACIATOSIAASGASGASGASTOSTITACATCOTGACTOSCOATOCASTACASCASGAACCASSGAACCASSTACTOTGACGASGAC ATOSTOCCOSSTICOLACCASSGACCASSTACASOTACCOCTTOSTCT ACIATOSIAAGGASGASGASTACTOTTACATCOTGACTOSCOATOCASTACASCASGAACCASSGAACCASSTACTOSTACACCASGAACCASGAACCASSTACCASGACCASSTACTOSTACASCASGAC ATOSTOCCOSSTICOLACCASSGACCASSTACASOTACCOCTTOSTCT ACIATOSIAAGGASGASGASTOSTITACATCOTGACCOSGATOCGASTACCASGASACCASSTACTOSGACGASGACCASSTACGACGASGACCASSTACTOSGACGASGACGASTGSTACACCTOSTACCACCOSTICOLACCASSGACCASSTACCASGACCASSTACGACGASGACGASGACGASTGSTACACCTOSTACCCCCGASTACCASGACCASSTACCASGACCASSTACCASGACCASSTACCASGACGASGACGASTGSTACACCTOSTACCACCOSGACCASSTACCASGACCASSTACCASGACCASSTACGACGASGACGASGACGASGACGASGACGASGACGASGACGASGACGASTGSTACCACCTOSGACCASSTACGACGASGACGASGACGASGACGASGACGASGACGASGACGASTGSTACCACCTOSTACCACCOSGACCASSTACCASGACCASSTACCASGACGACGASGACGACGACGACGACGACGACGACGACGACGACGACGACG
ITASRG ITA4PG Turkey Rio Mexico MLA MIP ITA3 ITA2 ITA1 Bb1NX-C4 GP-C19	ATOSTOCCOGITECOSACCAASGACCAASTACAGOTACCOCTEOSTET ACGATOGAAAGGAGAGTOTTIACATCETOTACTOSCOGITECASTACAACCAGOGAACCAGOTACCISTECTOTACTOSCOGITECOSACCAGAGGACCAASTACAGOTACCOCTEOSTET ACGATOGAAAGGAGGAGTOTTIACATCETOTACTOSCOGITECOSACCAGGGAACCAGOTACCOTTOSTET ACGATOGAAAGGAGGAGATOTOTTIACATCCTOTACTOSCOGITECOSACCAGGGACCAASTACAGOTACCOTTOSTET ACGATOGAAAGGAGGAGATOTOTTIACATCCTOTACTOSCOGITECOSACCAGGGACCAASTACAGOTACCOTTOSTET ACGATOGAAAGGAGGAGATOTOTTIACATCCTOTACTOSCOGITECOSACCAGGGACCAASTACAGOTACCOTTOSTET ACGATOGAAAGGAGGAGATOTOTTIACATCCTOTACTOSCOGITECOSACCAGGGACCAASTACAGOTACCOTTOSTET ACGATOGAAAGGAGGAGGAGGAGGAGGAGGAGGAGGAGGACGATOTTIACATCCTOSCOGITECOSACCAGGGACCAASTACAGOTACCOTTOSTET ACGATOGAAAGGAGGAGGAGGAGGAGGAGGAGGAGGAGGAGGAGG
ITASRG ITASRG Turkey Rio Mexico MIA MIP ITA3 ITA2 ITA3 ITA2 ITA1 BD1NW-C4	A TOSTOCCOSTITOCIA-CLABASSIACCASTACASOTACCOCTTOSTCT ACIATOSIAASSASAATOTOTTACATCOTTACTOSCOATOCASTACAACCASOSIAACCASOTACCOCTOSTCT ACIATOSIAASSASAATOTOTTACATCOTTACTOSCOATOCASTACAACCASOSIAACCASOTACCATOTISTCASOASSAC A TOSTOCCOSTITOCIA-CLABASSIACCASSTACASOTACCOCTTOSTCT ACIATOSIAASSASAATOTOTTACATCOTTOCTOSCOATOCASTACAACCASOSIAACCASOTACCATOSIACSASSAC A TOSTOCCOSTITOCIA-CLABASSIACCASSTACASOTACCOCTTOSTCT ACIATOSIAASSASAATOTOTTACATCOTTOCATCOSCOATOCASTACAACCASOSIAACCASOTACTOTGIACASOSAC A TOSTOCCOSTITOCIA-CLABASSACCASSTACASOTACCOCTTOSTCT ACIATOSIAASSASAATOTOTTACATCOTTOCATCOSCOATOCASTACAACCASOSIAACCASOSIACCASOTACTOTGIACCASOTAC
ITASRG ITA4PG Turkey Rio Mexico MLA MIP ITA3 ITA2 ITA1 Bb1NX-C4 GP-C19	A TOSTOCCOST TO CALCLARGE ACCASTA CAGOTA COCCTTOST CT ACCATOGAA GARGE AGATOST TACATOCTOTACTO COCCATOCASTA CALCCASOGAA CCASOTA CTOSTA COMO AGAS A TOSTOCCOST TO CALCLARGE ACCASTA CAGOTA COCCTTOST CT ACCATOGAA GARGE AGATOST TACATOCTOTACTO COCCATOCASTA CALCCASOGAA CCASOTA CTOSTO COLS SAG A TOSTOCCOST TO COLCARGE ACCASTA CAGOTA COCCTTOST CT ACCATOGAA GARGE AGATOST TACATOCTOTACTO COCCATOCASTA CALCCASOGAA CCASOTA CTOSTO COLS SAG A TOSTOCCOST TO COLCARGE ACCASTA CAGOTA COCCTTOST CT ACCATOGAA GARGE AGATOST TACATOCTOTACTO COCCATOCASTA CALCCASOGAA CCASOTA CTOSTO CAGOTACTO A TOSTOCCOST TO COLCARGE ACCASTA CAGOTA COCCTTOST CT ACCATOGAA GARGE AGATOST TACATOCTO TACCTO COCCATOCASTA CALCCASOGAA CCASOTA CTOSTO CAGOTACTO A TOSTOCCOST TO COLCARGE ACCASTA CAGOTA COCCTTOST CT ACCATOGAA GARGE AGATOST TACATOCTO TACCTO COCCATOCASTA CALCCASOGAA CCASOTA CTOSTO CAGOTACTO A TOSTOCCOST TO COLCARGE ACCASTA CAGOTA COCCTTOST CT ACCATOGAA GARGE AGATOST TACATOCTO TACCTO COCCATOCASTA CALCCASOGAA CCASOTA CTOSTO CAGOTACTO A TOSTOCCOST TO COLCARGE ACCASOTA COCCTTOST CT ACCATOGAA GARGE AGATOST TACATOCTO TACCTO COCATOCASTA CALCCASOGAA CCASOTA CTOSTO CAGOTACCASOTA COCCTOST CT ACCATOGA AGASGA GARGE AGATOST TACATOCTO TACCTO COCATOCASTA CALCCASOGAA CCASOTA CTOSTO CAGOTACTO A TOSTOCCOST TO COLCARGE AGAGA CALSTA CAGOTACCOCTTOST CT ACCATOGAA GARGE AGATOST TACATOCTO TACCTO COCATOCASTA CALCOSGAA CCASOTA CTOSTO CAGOTACTOST CACTOST CACCASOTA COCCTTOST CT ACCATOGA AGASGA GARGE AGATOST TACATOCTO TACCTO COCATOCASTA CALSOSTA CCASOTACTO TO CAGOTACCASOTA CCASOTA CCASOTACTO TO CACTOST CACCASOTACCASOTACTO TO CACTOST COCACTOST COCACTOST CACATOST CACATOST CACATOST CACTOST CACATOST CACATOST CACTOST CACATOST CACTOST CACATOST CACATOST CACTOST CACATOST CACAT
ITASRG ITAAPG Turkey Rio Mexico Mila MIP ITA3 ITA2 ITA1 B013W-C4 GP-C19 B01GP-C7	ATOSTOCCOGITECCALCABAGACCAGITACAGOTACCCCTTOSITE ACCATOGIAAGAGAGAGTITGTIACATCOTOTACTOSOCGATOCAGITACAACCAGGAACCAGGIACCAGITACAGCAGGAACCAGGIACCAGITACAGCAGGIACCACGICCAGGIACCAGGIACCACGICCAGGIACCACCAGGIAACCAGGIAACCAGGIACCAGGIA
ITASBG ITASBG TURKey RLO MEXICO MIA MIP ITA3 ITA3 ITA3 ITA3 ITA3 ITA3 ED15M-C4 GP-C19 ED1GP-C7 GP-C1	A TOSTOCCOST TO CALCLARGE ACCASTA CABOTACCOCTTOST CTACIA TOGALA GABORA ATTOST TACAT COTTACTO COCOLA TOCASTA CALCCASSA ACCASSTA CALCOTTOST CASTA CABORA CABOSA CABOTACTOST ACCASTA CALCASSA CABOTACTOST ACCASTG ACCASSA CABOTACTOST ACCASTG ACCASTA CABOTACCASSA CABOTACTOST ACCASTG ACCASSA CABOTACTOST ACCASTG ACCASTA CABOTACTOST ACCASTG ACCASSA CABOTACTOST ACCASSA CABOTACTOST ACCASTG ACCASSA CABOTACT ACTOST ACCASTG ACCASSA CABOTACT ACTOST ACCASSA CABOTACT ACCASTG ACCASSA CABOTACT ACTOST ACCASTG ACCASTA ACCASSA ACCASSA CABOTACT ACTOST ACCASTG ACCASSA ACCASSA ACCASSA ACCASSA ACCASSA ACCASSA CABOTACT ACTO
ITASRG ITASRG ITASRG Verkey Rio MMA MIP ITA3 ITA3 ITA3 ITA3 ITA3 ITA3 ITA3 ITA3	A TOSTOCCOST TO CALCAASSA CLASTACASSTACCOST TOST TA COATOSIAA GASGAGATOTTIACATOCTICACTOSOCIATOCASTACAACCASSAA CLASSTACTS TAGAGASASA A TOSTOCCOST TO COACCAASSACCAASTACASSTACCOST TOST TA COATOSIAA GASGAGATOTTIACATOCTISTACTOSOCIATOCASTACAACCASSAA CLASSTACTOSTIACASO A TOSTOCCOST TO COACCAASSACCAASTACASSTACCOST TOST TA COATOSIAA GASGAGATOTTIACATOCTISTACTOSOCIATOCASTACAACCASSGAA CLASSTACTISTIACASO A TOSTOCCOST TO COACCAASSGACCAASTACASSTACCOST TO TOT TA COATOSIAAA GASGAGATOTTIACATOCTISTACTOSOCIATOCASTACAACCASSGAA CLASSTACTISTIACASO A TOSTOCCOST TO COACCAASSGACCAASTACASSTACASSTACCOST TO STCA COATOSIAA GASGAGATISTI TA CATOCTISTACTOSOCIATI SCASTACAACCASSGAA CLASSTACTISTIACATOCTI A TOSTOCCOSST TO COACCAASSGACCAASTACASSTACCOST TO TOT TO COATOSIAAS GASGAGATISTI TA CATOCTISTACTOSOCIATI SCASTACAACCASSGAA CLASSTACTISTIGACIASAS A TOSTOCCOSST TO COACCASASGACCAASTACASSTACASSTACCOCTTOST CTACGATOSIAAGASGAGATISTI TA CATOCTISTACTOSOCIATI SCASTACAACCASSGAA CLASSTACTISTIGACIASASGAC A TOSTOCCOSST TO COACCASASGACCAASTACASSTACASSTACCOCTTOST CTACGATOSIAA GASGAGATISTI TA CATOCTISTACTOCOCIATI SCASTACAACCASSGAA CLASSTACTISTIGACGASSACCA A TOSTOCCOSST TO COACCASASGACCAASTACASSTACASSTACCOCTTOST CTACGATOSIAAGAGAGATISTI TA CATOCTISTACTOCOCIATI SCASTACCACCASSGAA CLASSTACTISTIGACGASSACCA A TOSTOCCOSST TO COACCASASGACCAASTACASSTACASSTACCOCTTOST CTACGATOSIAGAGAGAGATISTI TA CATOCTISTACATOCTISTICATOCCOCIATI SCASTACCACSASGAA CLASSTACTISTIGAGASSACCA A TOSTOCCOSST TO COACCASASGACCAASTACASSTACASSTACCOCTTOST CTACGATISTICATOCTISTACATOCTISTICATOCOCIATI SCASTACCASSGAA CLASSTACTISTIGAGASSACCA A TOSTOCCOSST TO COACCASASGACCAASTACASSTACCOCTTOST CTACGATOSI TA CATOCTISTACATOCTISTICATOCCOGATISCASTACCACSSGAA CLASSTACTISTIGAGASSACCA A TOSTOCCOSST TO COACCASASGACCAASTACASSTACCOCTTOST CTACGATOSI TACATOCTISTACATOCTISTICATOCCOGATISCASTACAACCASSGAA CLASSTACTISTIGACGASSACCA A TOSTOCCOSST TO COACCASASGACCAASTACASSTACCOCTTOST CTACGATOSI TACATOCTISTACATOCTISTICATOCCOGATISCASTACAACCASSGAA CLASSTACTISTIGACGASSACCA A TOSTOCCOSST TO COACCASASSACCAASTACASSTACCOCTTOST CTACGASGASGASGA
ITASSG ITA4PG Turkey Rio MMA NIP ITA3 ITA3 ITA3 ITA3 ITA3 ITA3 BD15N-C4 GP-C19 BD15N-C7 GP-C1 F3-109 F3-146	A TOSTOCCOST TO CALCLARGE ALCANST ALGO TACCOCTTOST CLACIA TOGALA GAGA GAGA TOTTI ALTOCT TOTAL TO COLOR TOLATACLACCASSA ALCASSA
ITASBG ITA4PG Urrkey RLO MMA MIP ITA3 ITA3 ITA3 ITA3 ITA3 ITA3 ITA3 ED15M-C4 GP-C19 ED16M-C7 GP-C19 ED16P-C7 GP-C19 F3-189 F3-146 WC-10284 WC-10284	A TOSTOCCOST TO CALCLARGE ALCANSTACABOTACCOCTTOST CTACIA TOSTA ACABORA ATTOST TACAT COTTO CALCOLOR TOCASTACAACCASOSA ACCASOSTACT STORA CALOSTACT TO CALCANSE ACCASOSTACT CONTRACT ACCAST TO CALCANSE ACCASOSTACT CONTRACT ACCAST TO CALCANSE ACCASOSTACT CONTRACT ACCAST TO CALCANSE ACCASOSTACT TO CALCANSE ACCASOST ACCASOST
ITASBG ITA4PG TUrkey Rio Mexico MiA MIP ITA3 ITA3 ITA3 ITA3 BD1MM-C4 GP-C19 BD1GP-C7 GP-C1 F5-189 F5-146 WC-90284 WC-951	A TOST OCCOSS TT COLACCABASIANCE AST ACASO TACCOCCTOST CT. ACIA TOSIA ACIASO AST ACTO TT ACT TOT TACT TO TOCOLA TOCAT TO ACIA CASO AL CONSTANCE TO TACA TO STANDARD AST ACIAS TACASO TACCOCCT TO ST C. ACIA TOSIA ACIASO AST ACIASO TACCOCT TO ST C. ACIA TOSIA ACIASO ASIA TO TA TACAT TO TO ACIA TOSIA ACIASO AL CASO TACT TO TACAT TO STANDARD AST ACIASO TACCOCT TO ST C. ACIA TOSIA ACIASO ASIA TO TA TACAT TO STANDARD ACIASO TACCOCT TO ST C. ACIA TOSIA ACIASO ASIA TO TA TACAT TO STANDARD ACIASO TACCATOSIA ACIASO TACCATOSIA ACIASO TACCATOSICA ACIASO AL CASO TACT TO STANDARD ACIASO TACCATOSIA ACIASO TACT TO STANDARD ACIASO TACT TO STANDARCASI ASIA CLASO TACT TO STANDARD ACIASO TACT TO STANDARCASI TACT TO STANDARCASI TO STANDARCASI TACATOSI TO STANDARCASI TACT TO STANDARCASI TACAT TO STANDARCASI TACT TO STANDAR
ITASBG ITA4PG Urrkey RLO MMA MIP ITA3 ITA3 ITA3 ITA3 ITA3 ITA3 ITA3 ED15M-C4 GP-C19 ED16M-C7 GP-C19 ED16P-C7 GP-C19 F3-189 F3-146 WC-10284 WC-10284	A TOSTOCCOST TO CALCLARAGE ACCASTA CABOTACCOCTTOST CTACIA TOGALA GRAGE AGOTA TOTTA CATO CTATACTOCOCA TOCASTA CALCLAGORA ACCASOTA CONSTRUCTA CON
ITASBG ITA4PG TUrkey Rio Mexico MIA MIP ITA3 ITA3 ITA3 ITA3 ITA3 ED107-C7 GP-C19 ED107-C7 GP-C19 ED107-C7 GP-C19 F5-146 WC-851 NP-C1244 WC-851 NP-C11	A TOSTOCCOST TO CALCAMBONICAMSTACABOTACCOCTTOST CTACIA TOSTA ACAGONA TOST TACATCOTTACATCOSCON TOCATACAACCASOBA E CANSTACTS TONA COMPANY A TOSTOCCOST TO CALCAMBONICAMSTACABOTACCOCTTOST CTACIA TOSTA AGANA SIA SON SATURDITACATCO TOTACTOSCON TOCATACAACCASOBA E CANSTACTS TONA COMPANY A TOSTOCCOST TO CALCAMBONICAMSTACABOTACCOCTTOST CTACIA TOSTA AGANO SAN STOTT TACATCO TOTACTOSCON TOCATACAACCASOBA E CANSTACTS TONA COMPANY A TOSTOCCOST TO COACCAMBONICAMSTACABOTACCOCTTOST CTACIA TOSTA AGANO SAN STOTT TACATCO TOTACTOSCON TO CALCASOBA E CANSTACTS TONA COMPANY A TOSTOCCOST TO COACCAMBONICAMSTACABOTACCOCTTOST CTACIA TOSTA AGANO SAN STOTT TACATCO TOTACTOSCON TO CALCASOBA E CANSTACTS TONA COMPONIC A TOSTOCCOST TO COACCAMBONICAMSTACABOTACCOCTTOST CTACIA TOSTA DAGANO SAN STOTT TACATCO TOTACTOSCON TO CALCASOBA E CANSTACTS TONA COMPONIC A TOSTOCCOST TO COACCAMBONICAMSTACABOTACCOCTTOST CTACIA TOSTA DAGANO SAN STOTT TACATCO TOTACTOSCON TO CALCASOBA E CANSTACT STONA COMPONIC A TOSTOCCOST TO COACCAMBONICAMSTACABOTACCOCTTOST CTACIA TOSTA DAGANO SAN STOTT TACATCO TOTACTOS CON TO CALCASOBA E CANSTACT STONA COMPONIC A TOSTOCCOST TO COACCAMBONICAMSTACABOTACCOCTTOST CTACINO SAN BONA SON SAN STOTT TACATCO TOTACTOS CON TO CALCASOBAN E CANSTACT STONA COMPONIC A TOSTOCCOST TO COACCAMBONICAMSTACABOTACCOCTTOST CTACINO SAN BONA SON SAN STOTT TACATCO TOTACTOS CON TO CALCASOBAN E CANSTACT STONA COMPONIC A TOSTOCCOSST TO COACCAMBONICAMSTACABOTACCOCTTOST CTACINO SAN BONA SON SAN STOTT TACATCO TOTACTOS CON TO CALCASOBAN E CANSTACT STONA COMPONIC A TOSTOCCOSST TO COACCAMBONICAMSTACABOTACCOCTTOST CTACINO SAN BONA SON SAN STOTT TACATCO TOTACTOS CON TO CANTACAACCASOBAN E CANSTACT STONA COMPONIC A TOSTOCCOSST TO COACCAMBONICAMSTACABOTACCOCTTOST CTACINO SAN BONA SON SAN STOTT TACATCO TOTACTOS CON TO CANTACAACCASOBAN E CANSTACT STONA COMPONIC A TOSTOCCOSST TO COACCAMBONICAMSTACABOTACCOCTTOST CTACINO SAN BONA SON SAN STOTT TACATCO TOTACTOS CON TO CANTACAACCASOBAN E CANSTACT STONA COMPONIC A TOSTOCCOSST TO COACCAMBONICAMSTACABOTACCOCTTOST CTACINO SAN BONA SON
ITASBG ITA4PG UTAPFG Verkey Rio MMA UTA3 ITA3 ITA3 ITA3 ITA3 BD15M-C4 GP-C19 BD15M-C4 GP-C19 BD15M-C7 GP-C1 F5-146 MC-10284 MP-C11 MP-C11 MP-C1	A TOST OCCOST TO CALCLARAGE ACCAST ACAGOTACCOCTTOST TA COLTOGRADA GRAGE AGOTATOT TACAT COTT ACCORDANCE AGOTACCOST OLDAN CAGOGAL CLASS TACAGOTACCOCTTOST TA COLTOGRADA GRAGE AGOTATOT TACAT COTT ACCORDANCE AGOTACCOST OLDAN CLASS TACAGOTAL COCTOST TO COLTOGRADA GRAGE AGOTACTOT TACAT COTTACCT COSC DA TOCAST ACAAC CLASS TACAGOTAL COTTOST TO COLTOGRADA GRAGE AGOTACTOST DA COLT
ITASBG ITA4PG TUrkey Rio Mexico MIA MIP ITA3 ITA3 ITA3 ITA3 ITA3 ITA3 ITA3 ITA3	A TOST OCCOST TO CALCLARRIGALCARSTALAGOTA COCCTTOST CTACIA TOGALAGAGAGATOTTILA. TO CITATA COCCIA TOLASTALAACAGOGAA COLOGITA COLASTALAACAGOTA COCCTTOST CTACIA TOGALAGAGAGATOTTILA. TO CITATA COCCIA TOLASTALAACAGOGAA COLOGITACIA COLASTALAGAGAGATACIA TO CITALA COCTIA COLASTALAACAGOGAA COLOGITACIA COLASTALAGAGAGAGATACIA TO CITALA COCTIA COLASTALAGAGAGAGATACIA TO CITALA COCTIA COLASTALAGACAGOTAL COCTIGOT CACIA TOGALAGAGAGAGATOTTILA. TO CITALA COLASTALAGACAGOGAA COLASTALAGATACIA COLASTALAGAGAGAGAGAGAGAGATOTTILA. TO CITALA COLASTALAGACAGOGAA COLASTALAGACAGOGAA COLASTALAGACAGOTAL COCTIGOT CACIA TOGALAGAGAGAGAGAGAGAGATOTTILA. TO CITALAGT COCCANTICA CACIASGAA COLASGAA COLASTALAGACAGOGAA COLASTALAGACAGOGAA COLASTALAGACAGAGAGAGAGAGAGAGAGAGAGAGATOTTILA. TO CITALAGT COCCANTICA CAGOGAA COLASGAA COLAGOTACIA COLASGAACAGAGAGAGAGAGAGAGAGAGAGAGAGAGAGAGAG
ITASBG ITA4PG UTXey Rio Mexico MiA MIP ITA3 ITA3 ITA3 ITA3 ITA3 BD10P-C7 GP-C1 9B10P-C7 GP-C1 9B10P-C7 GP-C1 95-146 WC-10284 WC-01284 WP-C11 MP-C1 DUVA KC3-C24	A TOSTOCCOGITICOLACIABIGAICA STA CAGOTACCOCTTOSTIC A COLTOGRADA GUAGA DITOTTA CATOCITAL TOCOCOLI SOCATACIAA CLAGGA A CLAGGTACISTICA COLOGIA DI CASTACIAA CLAGGIA CONSTRUCTASI DI CASTACIAA CLAGGIA

Multiple alignment of a conserved region of the *ama-1* gene of the strains of *Babesia bigemina* reported in the GenBank. The asterisks () show the conserved regions, the blank spaces show the nucleotide differences of the analyzed sequences, and the colored letters show that the nucleotides vary in one or two strains.

have the same level of detection, generating an amplification in a dilution of 1×10^{-8} (Figure 6).

Evaluation of the LAMP test with DNA from *B. bigemina* strains from different geographical areas


Technique validation was carried out using DNA samples of *B. bigemina* from different geographical areas, such as Chiapas, Colima, Jalisco, Nayarit, Queretaro, Colombia, and the Argentina strains, S1A, S2P and S3P. The results show that the characteristic amplicons indicate the amplification of the target sequence in all analyzed strains. Additionally, the colorimetric result, consistent with the electrophoretic results, showed the color change from violet to blue, indicating a positive result (Figure 7).

Discussion

Bovine babesiosis is a disease of cattle whose timely diagnosis is an important tool for controlling and preventing the spread of the disease. The use of effective techniques is aimed at covering five general purposes: (1) identification of blood parasite species and strains involved in the disease, (2) determination of the distribution and evaluation of disease risk in a herd at the


Primer	Size (bp)	Sequence (5'-3')
F3	20	TGTGATCAAGCACACTGGAG
B3	20	CCATCGTAGACGAAGGGGTA
FIP	39	ACCCCACTTACGGAGGTCGGGCTCTTACAAACGCAGGCA
BIP	43	GAAACGCCGTTACTAACTGCGCTGTACTTGGTCCTCTGGTCG
LF	21	GTAGGTGCCATTAAACAGCGG
LR	17	CAGCAACATCGTGCCCG





regional or national level, (4) certification of the state of infection of the animals by commercial requirements or in eradication programs, (5) identification or confirmation of the cause of outbreaks and deaths, and (6) identification of specific arthropods as vectors or of the pathogens transmitted by them (25).

The LAMP test has been considered a novel detection technology and is widely used because of its competitive advantages. First, when compared with traditional detection methods, the specificity of LAMP is high due to the quantity of oligonucleotides that are required. The six (or eight) oligonucleotides are specifically designed to recognize six (or eight) different regions on the target DNA. If one of the oligonucleotides does not hybridize correctly, the reaction does not occur, and the fact that all oligonucleotides are required to hybridize eliminates non-specific amplification. Second, LAMP




has been shown to be as sensitive as nPCR and is capable of detecting 10 or even fewer copies of template DNA (14, 26–29). That is, the test can be used to detect the pathogen's DNA in early stages of infection in animals with no apparent clinical signs or in carrier animals. In this work, the molecular detection of *B. bigemina* was evaluated using the LAMP test by detecting and amplifying a conserved region of the *ama-1* gene. The *ama-1* gene of *B. bigemina*, initially characterized by Torina et al. (17) showed a high degree of similarity between the nucleotide sequences available in the Genbank (99.98%). Subsequent research has used this gene as a detection target gene in epidemiological studies to detect *B. bigemina* (30, 31), suggesting high specificity in molecular tests in addition to a high degree of conservation between isolates. Because species-specific detection tools are of great relevance, especially

infections with closely related species of the parasite exist, detection techniques that use ama-1 sequences as detection targets are an excellent alternative (30). The LAMP test was found to be specific for the detection of B. bigemina as shown by bioinformatic methods and experimental verification since no amplification or color changes were observed, indicating a positive reaction to the DNA from other species. In the specificity assay, we did not include other Babesia or Theileria species because we did not have them at the time and more importantly, in the primer design analyses, no similarity was found between the four primers designed with any sequence published in the GenBank from other Babesia or Theileria species. However, biological confirmation must be done in the future including DNA samples from Theileria and Babesia species that infect cattle. To date, several studies carried out on the sequence of the ama-1 gene characterize it as a conserved gene between isolates and strains and divergent from other related species, which makes it an excellent candidate for the development of techniques that use sequences of this gene as a molecular detection target (30, 32). Additionally, this study shows that the LAMP test can detect B. bigemina strains from different regions of Mexico (Chiapas, Colima, Jalisco, Nayarit, and Queretaro), and from two South American countries, namely Argentina (three strains) and Colombia. In previous reports concerning the LAMP technique, the authors did not evaluate the detection of different strains for which the present work constitutes the first LAMP test that demonstrates the detection of B. bigemina from different geographical regions. These results suggest that the LAMP test for the detection of B. bigemina can be implemented in different parts of the world.

One of the essential elements to develop a molecular, sensitive, and specific detection test using the LAMP technique is the design of oligonucleotides (26). Previous studies have shown that the design determines the success of the technique, paying attention to the formation of dimers, and although primer purification is not a strict requirement, it is a recommended step to increase the specificity and reproducibility of the test (33). As a note of caution, tests such as nested PCR, qPCR, or LAMP, which are ultra-sensitive, can have contamination issues if proper handling, conditions, and/ or facilities are lacking. For example, sample reception area, nucleic acid extraction area, and analysis areas are required. Ideally, these areas should be physically separated.

Within the standardized factors, a differential point of the present investigation was the elimination of the denaturation step (11). Although this factor was not compared with other methods, the reaction did not need this step to be carried out. On the other hand, the isothermal amplification temperature is decisive for performing the technique without using specialized equipment, such as a thermocycler. The results indicate that the amplification proceeded at all evaluated temperatures. This process provides a competitive advantage by offering flexibility



in this factor since it allows thermal variation between the equipment used for the test.

Among the advantages of the test is the diagnostic specificity. The results from our study show that the detection method proposed in this study is highly specific for the diagnosis of B. bigemina, which gives it an added advantage over the previously developed nPCR test (24) since the latter, in addition to B. bigemina, also amplifies B. ovata DNA. Bioinformatic analysis of the primers designed in this work did not identify the sequence of any gene in B. ovata including ama-1. This study is the first one that uses the conserved region of the ama-1 gene for the molecular diagnosis using the LAMP technique of bovine babesiosis caused by B. bigemina. Additionally, it is the first study to consider that within its bioinformatic analysis, strains of Babesia bigemina from all over the world as reported in the GenBank can be analyzed. Likewise, strains from different parts of the world were experimentally evaluated, which suggests that this technique could have a worldwide application.

On the other hand, the diagnostic sensitivity showed that the LAMP technique has the same detection limits as nPCR. Other studies concerning the LAMP technique show that the technique is 10 times more sensitive than nPCR and 100 times more sensitive than endpoint PCR (10, 15, 16). In these studies, different target genes and different complementary methodology test were used. Additionally, the serial dilutions of the DNA were performed in different ways in each of the aforementioned studies, accurately quantifying the DNA present in each dilution. In other studies, they diluted infected

07






FIGURE 5

Evaluation of the diagnostic sensitivity of the LAMP technique, for the detection of the *ama-1* gene of *B. bigemina*. (A) Electrophoresis in 1.2% agarose gel, stained with red gel. (B) Colorimetric test with hydroxynaphthol blue. Lanes: M) Molecular size marker in bp; (1) Negative control, (2) 2% parasitaemia, (3) 1% parasitaemia, (4) 1×10^{-1} , (5) 1×10^{-2} , (6) 1×10^{-3} , (7) 1×10^{-4} , (8) 1×10^{-5} , (9) 1×10^{-6} , (10) 1×10^{-7} , (11) 1×10^{-8} , and (12) 1×10^{-9} .



erythrocytes with uninfected erythrocytes to reach a known concentration of parasites (10, 15, 16). In this study, sensitivity was shown when determining the percentage of parasitemia, obtaining a sensitivity of 0.00000001% of parasitemia thus demonstrating a high sensitivity. In addition, compared to the nPCR technique, which is considered the gold standard, the advantages offered by LAMP are low cost (almost half the cost of nPCR), reduced time (2 h for LAMP vs. >5 h

for nPCR), accessibility of the technique, high sensitivity, and specificity. Finally, one of the most important qualities of LAMP as shown in the present study is that it makes it possible to interpret the results visually through colorimetry. The use of hydroxynaphthol blue (HNB) gives it superiority because in addition to being economical, it is not necessary to open the reaction tube after the incubation time, a step that could lead to contamination by the amplification products



(13, 26, 33). Furthermore, as mentioned by Wastling et al. (34), who compared various visual detection methods, HNB is the best visual detection method due to its ease of application, cost, and sensitivity. In conclusion, the LAMP technique for the diagnosis of *B. bigemina* when compared with the nPCR technique, which is considered the gold standard, offers advantages, such as low cost, reduced time, accessibility of the technique, and both high sensitivity and specificity. Comparison with other gold standard techniques and further steps in the validation process using field samples are underway, including the detection of *B. bigemina* in infected ticks.

Data availability statement

The original contributions presented in the study are included in the article/supplementary materials, further inquiries can be directed to the corresponding author/s.

References

1. Babes V. Sur l'hemoglobinurie bacterienne du boeuf. CR Acad Sci. (1888) 107:692-4.

2. Bock R, Jackson L, Vos AD, Jorgensen W. Babesiosis of cattle. *Parasitology*. (2004) 129:S247-69. doi: 10.1017/S003118200400 5190

Ethics statement

The animal study was reviewed and approved by Bioethics Committee of the College of Natural Sciences, Autonomous University of Queretaro, Mexico (Approval Numbers; FCN/2011-0221 and 71FCN2016).

Author contributions

AL-Z conducted the experiments, analyzed the results and wrote the initial manuscript. BC-G participated in the experimental design and wrote the manuscript. SW contributed with biological DNA samples, reagents, and protocol design. SC helped design the primers and standardize the test. MT helped design the primers and optimize the test. JM conceptualized the idea, obtained the funds, and edited the final version of the manuscript. All authors revised and approved the final version of the manuscript.

Funding

This research was funded by CONACyT-Mexico Problemas Nacionales 2017 and by the USA Agricultural Research Service 59-2090-1-001-F.

Conflict of interest

The authors declare that the research was conducted in the absence of any commercial or financial relationships that could be construed as a potential conflict of interest.

Publisher's note

All claims expressed in this article are solely those of the authors and do not necessarily represent those of their affiliated organizations, or those of the publisher, the editors and the reviewers. Any product that may be evaluated in this article, or claim that may be made by its manufacturer, is not guaranteed or endorsed by the publisher.

^{3.} Smith T, Kilborne FL. Investigations Into the Nature, Causation, and Prevention of Texas or Southern Cattle Fever. (1893). Available online at: http://agris.fao.org/ agris-search/search.do?recordID=US201300655646 (accessed February 11, 2022). doi: 10.5962/bhl.title.124068

^{4.} Yu L, Terkawi MA, Cruz-Flores MJ, Claveria FG, Aboge GO, Yamagishi J, et al. Epidemiological Survey of *Babesia bovis* and *Babesia bigemina* Infections of Cattle in Philippines. *J Vet Med Sci.* (2013) 75:995–8. doi: 10.1292/jvms.12-0425

^{5.} Schnittger L, Rodriguez AE, Florin-Christensen M, Morrison DA. *Babesia*: a world emerging. *Infect Genet*. (2012) E12:1788–809. doi: 10.1016/j.meegid.2012.07.004

^{6.} Smith RD, Evans DE, Martins JR, Ceresér VH, Correa BL, Petraccia C, et al. Babesiosis (*Babesia bovis*) stability in unstable environments. *Ann N Y Acad Sci.* (2000) 916:510–20. doi: 10.1111/j.1749-6632.2000.tb05330.x

7. Suarez CE, Noh S. Emerging perspectives in the research of bovine babesiosis and anaplasmosis. *Vet Parasitol.* (2011) 180:109–25. doi: 10.1016/j.vetpar.2011.05.032

8. Zintl A, Gray JS, Skerrett HE, Mulcahy G. Possible mechanisms underlying age-related resistance to bovine babesiosis. *Parasite Immunol.* (2005) 27:115–20. doi: 10.1111/j.1365-3024.2005.00748.x

9. Criado-Fornelio A, Martinez-Marcos A, Buling-Saraña A, Barba-Carretero JC. Molecular studies on *Babesia*, *Theileria* and *Hepatozoon* in southern Europe. *Part I Epizootiological aspects Vet Parasitol.* (2003) 113:189–201. doi: 10.1016/S0304-4017(03)00078-5

10. Iseki H, Alhassan A, Ohta N, Thekisoe OMM, Yokoyama N, Inoue N, et al. Development of a multiplex loop-mediated isothermal amplification (mLAMP) method for the simultaneous detection of bovine *Babesia* parasites. *J Microbiol Methods*. (2007) 71:281–7. doi: 10.1016/j.mimet.2007.09.019

11. Notomi T, Okayama H, Masubuchi H, Yonekawa T, Watanabe K, Amino N, et al. Loop-mediated isothermal amplification of DNA. *Nucleic Acids Res.* (2000) 28:E63. doi: 10.1093/nar/28.12.e63

12. Nagamine K, Hase T, Notomi T. Accelerated reaction by loop-mediated isothermal amplification using loop primers. *Mol Cell Probes.* (2002) 16:223–9. doi: 10.1006/mcpr.2002.0415

13. Goto M, Hond, E, Ogura A, Nomoto A, Hanaki KI. Colorimetric detection of loop-mediated isothermal amplification reaction by using hydroxy naphthol blue. *BioTechniques*. (2009) 46:167–72. doi: 10.2144/000113072

14. Parida M, Sannarangaiah S, Dash PK, Rao PV, Morita K. Loop mediated isothermal amplification (LAMP): a new generation of innovative gene amplification technique; perspectives in clinical diagnosis of infectious diseases. *Rev Med Virol.* (2008) 18:407–21. doi: 10.1002/rmv.593

15. Liu A, Guan G, Du P, Gou H, Liu Z, Liu J, et al. Loop-mediated isothermal amplification (LAMP) method based on two species-specific primer sets for the rapid identification of Chinese *Babesia bovis* and B. *bigemina. Parasitol Int.* (2012) 61:658–63. doi: 10.1016/j.parint.2012.07.004

16. Yang Y, Li Q, Wang S, Chen X, Du A. Rapid and sensitive detection of *Babesia bovis* and *Babesia bigemina* by loop-mediated isothermal amplification combined with a lateral flow dipstick. *Vet Parasitol.* (2016) 219:71–6. doi: 10.1016/j.vetpar.2016.02.004

17. Torina A, Agnone A, Sireci G, Mosqueda JJ, Blanda V, Albanese I, et al. Characterization of the apical membrane antigen-1 in Italian strains of *Babesia bigemina*. *Transbound Emerg Dis.* (2010) 57:52-6. doi: 10.1111/j.1865-1682.2010.01118.x

18. Moreau E, Bonsergent C, Al Dybiat I, Gonzalez LM, Lobo CA, Montero E, et al. *Babesia divergens* apical membrane antigen-1 (BdAMA-1): a poorly polymorphic protein that induces a weak and late immune response. *Exp Parasitol.* (2015) 155:40–5. doi: 10.1016/j.exppara.2015.04.024

19. Rittipornlertrak A, Nambooppha B, Simking P, Punyapornwithaya V, Tiwananthagorn S, Jittapalapong S, et al. Low levels of genetic diversity associated with evidence of negative selection on the *Babesia bovis* apical membrane antigen 1 from parasite populations in Thailand. *Infect Genet.* (2017) E54:447–54. doi: 10.1016/j.meegid.2017.08.009

20. Centro Nacional para la Información Biotecnológica. NCBI. GenBank. (2016). Available online at: https://www.ncbi.nlm.nih.gov/genbank/ (accessed September 22, 2022).

21. PrimerExplorer. Eiken Chemical Co., Ltd. Available online at https://primerexplorer.jp/e/index.html (accessed May 29, 2022).

22. Eiken Chemical, Co., Ltd. A Guide to LAMP Primer Designing. (2019). Available online at: https://primerexplorer.jp/e/v5_manual/pdf/ PrimerExplorerV5_Manual_1.pdf. (accessed May 29, 2022).

23. Organización Mundial, d. e. Sanidad Animal. OIE. Sistema Mundial de Información zoosanitara. WAHIS-Wild Interface. (2015). Available online at: https://www.woah.org/fileadmin/Home/eng/Health_standards/tahm/3.04.02_BABESIOSIS.pdf (accession February 11, 2022).

24. Figueroa JV, Chieves LP, Johnson GS, Buening GM. Multiplex polymerase chain reaction-based assay for the detection of *Babesia bigemina*, *Babesia bovis* and *Anaplasma marginale* DNA in bovine blood. *Vet Parasitol.* (1993) 50:69–81. doi: 10.1016/0304-4017(93)90008-B

25. Böse R, Jorgensen WK, Dalgliesh RJ, Friedhoff KT, Vos AJ. Current state and future trends in the diagnosis of babesiosis. *Vet Parasitol.* (1995) 57:61– 74. doi: 10.1016/0304-4017(94)03111-9

26. Tomita N, Mori Y, Kanda H, Notomi T. Loop-mediated isothermal amplification (LAMP) of gene sequences and simple visual detection of products. *Nat Protoc.* (2008) 3:877–82. doi: 10.1038/nprot.2008.57

27. Parida M, Posadas G, Inoue S, Hasebe F, Morita K. Real-time reverse transcription loop-mediated isothermal amplification for rapid detection of West Nile virus. *J Clin Microbiol.* (2004) 42:257–63. doi: 10.1128/JCM.42.1.257-263. 2004

28. Parida M, Horioke K, Ishida H, Dash PK, Saxena P, Jana AM, et al. Rapid detection and differentiation of dengue virus serotypes by a real-time reverse transcription-loop-mediated isothermal amplification assay. *J Clin Microbiol.* (2005) 43:2895–903. doi: 10.1128/JCM.43.6.2895-2903.2005

29. Hong TC, Mai QL, Cuong DV, Parida M, Minekawa H, Notomi T, et al. Development and evaluation of a novel loop-mediated isothermal amplification method for rapid detection of severe acute respiratory syndrome coronavirus. *J Clin Microbiol.* (2004) 42:1956–61. doi: 10.1128/JCM.42.5.1956-1961. 2004

30. Sivakumar, T., Altangerel, K., Battsetseg, B., Battur, B., Aboulaila, M., Munkhjargal, T., et al. Genetic detection of *Babesia bigemina* from Mongolian cattle using apical membrane antigen-1 gene-based PCR assay. *Vet Parasitol.* (2012) 8:187. doi: 10.1016/j.vetpar.2012.01.008

31. Elsify A, Sivakumar T, Nayel M, Salama A, Elkhtam A, Rizk M, et al. An epidemiological survey of bovine *Babesia* and *Theileria* parasites in cattle, buffaloes, and sheep in Egypt. *Parasitol Int.* (2015) 64:79–85. doi: 10.1016/j.parint.2014.10.002

32. Ganguly A, Maharana BR, Ganguly I. Pentaplex PCR assay for rapid differential detection of *Babesia bigemina*, *Theileria annulata*, *Anaplasma marginale* and *Trypanosoma evansi* in cattle. *Biologicals*. (2020) 63:81–8. doi: 10.1016/j.biologicals.2019.10.011

33. Tanner NA, Zhang Y, Evans TC. Simultaneous multiple target detection in real-time loop-mediated isothermal amplification. *Biotechniques*. (2012) 53:81–9. doi: 10.2144/0000113902

34. Wastling SL, Picozzi K, Kakembo ASL, Welburn SC. LAMP for human African trypanosomiasis: a comparative study of detection formats. *PLoS Negl Trop Dis.* (2010) 4:e865. doi: 10.1371/journal.pntd.000 0865